**REF 50-7050**

## **PleurX™ Pleural Catheter Mini Kit**

- en** Pleural Catheter Mini Kit  
For Pleural Placement Only
- fr** Mini kit cathéter pleural  
Placement pleural uniquement
- de** Pleurakatheter-Mini-Kit  
Nur zur pleuralen Anwendung
- it** Mini kit per catetere pleurico  
Esclusivamente per  
posizionamento pleurico
- es** Mini kit del catéter pleural  
Solo para colocación pleural
- pt** Mini kit de cateter pleural  
Apenas para aplicação pleural
- nl** Minikit voor pleurale katheter  
Alleen voor pleurale plaatsing
- sv** Pleurakateter, miniset  
Endast för pleural placering
- fi** Keuhkopussikatetriminiseti  
vain keuhkopussiin asettamiseen
- no** Pleurakateter – minisett  
Bare for pleural plassering
- da** Pleurakateter, minisæt  
Kun til anlæggelse i pleura
- el** Μίνι κιτ καθετήρα υπεζωκοτικής κοιλότητας  
Για τοποθέτηση στην υπεζωκοτική κοιλότητα μόνο
- tr** Plevral Kateter Mini Kit  
Yalnız Plevral Boşluğa Yerleştirilmeye Yöneliktir

Not made with  
natural rubber  
latex.

**STERILE EO**

For Single Use Only

Do Not Resterilize

Do not use if  
package is damaged.

**DEHP**  
Contains phthalates

**CareFusion**  
75 North Fairway Drive  
Vernon Hills, IL 60061 USA

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Sponsored by  
**CareFusion Australia 316 Pty. Ltd.**  
**CareFusion New Zealand 313 Ltd.**

361-34701 • 2017-05

**USA Rx Only**

© 2010, 2011, 2012, 2013, 2014, 2015, 2017,  
CareFusion Corporation or one of its affiliates.  
All rights reserved.

For product inquiries or technical assistance:  
pleurx-info@bd.com www.bd.com/pleurx

PleurX, CareFusion and the CareFusion logo  
are trademarks or registered trademarks of  
CareFusion Corporation, or one of its affiliates.

BioPatch and Dermabond are trademarks  
of Johnson & Johnson Corporation.

DCN 60575-2  
Vernon Hills, IL  
Richard Cisneroz  
05-10-17  
361-34701

Proofed by: \_\_\_\_\_ Date: \_\_\_\_\_  
Dimensions checked: \_\_\_\_\_ Copy checked: \_\_\_\_\_









Figure, Abbildung, Figura, Afbeelding, Figur, Kuva, Εικόνα, Şekil


DCN 60575-2  
Vernon Hills, IL  
Richard Cisneroz  
05-10-17  
361-34701

Proofed by: \_\_\_\_\_ Date: \_\_\_\_\_  
Dimensions checked: \_\_\_\_\_ Copy checked: \_\_\_\_\_


Figure, Abbildung, Figura, Afbeelding, Figur, Kuva, Εικόνα, Şekil


**Figure, Abbildung, Figura, Afbeelding, Figur, Kuva, Εικόνα, Şekil**


**(12)**


**(17)**

**(13)**


**(18)**


"CLIQUEZ" "KLIČK" "CLIC" "KLIK"  
 "CLIQUE" "KLIČK PÁ" "NAKSAHDUS"  
 "KLIK" "Κλικ"


**(14)**


**(19)**


**(15)**

**(20)**

**(16)**



**(21)**

## PleurX Pleural Catheter Procedure

### Product Information

The PleurX Pleural Catheter consists of a fenestrated silicone catheter with a valve mechanism and a polyester cuff. A barium sulfate stripe runs the entire length of the catheter. (Figure 1) The valve is designed to prevent the passage of air or fluid in either direction unless it is accessed with the specifically matched drainage line, access kit, or vacuum bottles provided by CareFusion. The PleurX Pleural Catheter is designed exclusively for use with the PleurX Vacuum Bottles, glass vacuum bottles, and the PleurX Lockable Drainage Line for connection to wall suction or portable suction.

### Indications for Use

The PleurX Pleural Catheter Mini Kit is indicated for intermittent, long term drainage of symptomatic, recurrent, pleural effusion, including malignant pleural effusion and other recurrent effusions that do not respond to medical management of the underlying disease. The devices are indicated for the palliation of dyspnea due to pleural effusion and providing pleurodesis (resolution of the pleural effusion). The PleurX Pleural Catheter is indicated for adults only.

The PleurX Lockable Drainage Line (Figure 2) is used to drain fluid using standard wall suction, water seal drainage system, vacuum bottle, or other appropriate method.

### Contraindications

Use of the PleurX Pleural Catheter Mini Kit is contraindicated in the following situations:

1. When there is a shift  $\geq 2$  cm in the mediastinum towards the ipsilateral side of the effusion.
2. When the pleural cavity is multi-loculated, and the drainage of a single loculation would not be expected to provide relief of dyspnea.
3. When there is a coagulopathy.
4. When the pleural cavity is infected.
5. When the effusion is known to be chylous.
6. When the patient is known or suspected to be allergic to materials contained in the device.

### Warnings

Do not put anything except the access tip of the PleurX Lockable Drainage Line, PleurX Catheter Access Kit, or PleurX Vacuum Bottles into the PleurX Catheter valve since any other device could damage the valve. A damaged valve may allow air into the body or let fluid leak out through the valve when not draining.

### Cautions

For single use only. Re-use may result in a non-functional product or contribute to cross contamination.

Sterile technique should be used when placing and draining the catheter.

Exercise care when placing the catheter to prevent it from coming into contact with surfaces such as drapes or towels. Silicone rubber is highly electrostatic and attracts airborne particles and surface contaminants.

The fenestrations must be entirely within the pleural space to avoid leakage into the tunnel tract. Take patient size, tunnel length, and the catheter length into account when placing the catheter.

Individuals with severe skin pathology overlying the favored implantation site (e.g. tumor invasion/psoriasis) should receive alternative treatment to relieve symptoms.

### Sterility

This product has been sterilized. It is for single use only and is not to be re-sterilized. Do not use if package is damaged. CareFusion will not be responsible for any product that is re-sterilized, nor accept for credit or exchange, any product that has been opened but not used.

### PleurX Catheter Tray Supplies

#### Placement Components

- 1 PleurX Pleural Catheter, 15.5 Fr
- 1 Guidewire Introducer with Needle
- 1 Syringe, 10 ml
- 1 J-Tip Guidewire
- 1 Peel-Away Introducer, 16 Fr
- 1 Tunneler
- 1 Needle Foam Stop
- 1 Blue Wrap

#### Drainage Components

- 1 Drainage Line with Lockable Access Tip
- 1 Needle, 17 G x 1" (2.5 cm)
- 1 5-in-1 Drainage Line Adapter
- 1 Valve Cap

#### Dressing Components


- 6 Gauze Pads, 4" x 4" (10.2 cm x 10.2 cm)
- 1 Foam Catheter Pad
- 1 Self-Adhesive Dressing

### General Guidelines

1. The procedure for pleural placement can be performed using local anesthetic and sedation. However, depending on patient needs, it may be performed using alternative approaches to anesthesia or sedation.

- Use of image guidance may aid in the accuracy and safety of catheter placement. Care should be taken to identify and avoid contact with vasculature near the guidewire insertion site.
- Guidewire insertion site selection should be based upon patient anatomy and presentation with consideration given to any possible adhesions or loculated pockets of fluid.
- Consideration should be given to the patient's ease of access in determining the location of the catheter exit site.
- The fenestrated end of the catheter may be cut shorter depending on an individual patient's anatomy. **Note:** Leave at least one fenestration present on the catheter.
- When using the drainage line, it may be used to drain at a suction level of -60 cm H<sub>2</sub>O or less for as long as it takes to drain 1,000 ml of fluid from the chest per day. The volume of fluid removed should be based on the individual patient's status and the risks for over draining, including re-expansion pulmonary edema in the chest.

### Suggested Placement Procedure

(3)

Proper medical and surgical procedures are the responsibility of the physician. The appropriateness of any procedure must be based upon the needs of the patient. **Figure (3)** illustrates the placement of the PleurX Pleural Catheter, as described in the following procedure.

- Position the patient appropriately to access the desired guidewire insertion site.
- Identify the appropriate intercostal space for guidewire placement. The guidewire is typically placed in the sixth or seventh intercostal space. Ultrasound can be used to confirm the guidewire insertion site.
- Identify the location of the catheter exit site, which is usually approximately 5 cm from the guidewire insertion site.
- Surgically prep the patient.

- Drape and anesthetize the planned insertion and tunneling sites.

**Caution:** Care must be taken when inserting the needle to avoid puncturing or lacerating the lung or liver.

- Insert the guidewire introducer with needle, attached to a syringe, through the desired intercostal space and just over the lower rib.
- Ensure free aspiration of pleural fluid, then remove the needle and syringe, leaving the guidewire introducer in place.

- Insert the guidewire through the introducer, advancing it well into the pleural cavity.

**Caution:** Do not allow the guidewire to inadvertently advance totally within the patient. Ensure that the guidewire exits the proximal end of the needle, dilator, or sheath prior to and during placement of each component.

- Remove the introducer, leaving the guidewire in place.

**Caution:** Damage to the guidewire may result if withdrawn through the needle.

- Make a 1 cm incision at the guidewire insertion site.

- Make a second 1-2 cm incision approximately 5 cm from the guidewire insertion site. This incision will be the catheter exit site. Consider the patient's ease of access in determining its location. **(Figure 3)**

**Note:** A smaller incision may provide better security of the catheter.

**Note:** Take care to ensure that the tunnel track has been anesthetized.

- The fenestrated end of the catheter may be cut shorter depending on an individual patient's anatomy. If desired, cut a portion of the fenestrated end of the catheter by using a scalpel to make a straight cut between fenestrations. Leave at least one (1) fenestration on the catheter.

- Attach the fenestrated end of the catheter onto the tunneler.

**Caution:** Exercise care when handling the catheter to prevent it from coming into contact with surfaces such as drapes or towels. Silicone rubber is highly electrostatic and attracts airborne particles and surface contaminates.

**Caution:** Use rubber-shod instruments when handling the catheter. Possible cuts or tears can occur if rubber-shod instruments are not used.

- Pass the tunneler **(A)** and catheter **(B)** subcutaneously from the second incision up to and out through the first incision at the guidewire insertion site. **(Figure 4)** Continue to draw the catheter through the tunnel until the polyester cuff lies inside the tunnel, about 1 cm **(C)** from the second incision. **(Figure 5)** Disconnect the tunneler from the catheter.

**Note:** If the cuff is advanced further into the tunnel, it can make later removal of the catheter difficult.

15. Thread the 16 Fr peel-away introducer sheath over the guidewire into the pleural cavity.
16. Remove the guidewire and dilator as a unit, leaving the 16 Fr peel-away introducer sheath in place.

**Caution:** Place a thumb over the end of the sheath as the dilator is removed to avoid air entering the pleural cavity. Care must be taken not to bend or kink the sheath. Damage to the sheath may prevent passage of the catheter.

17. Insert the fenestrated end of the catheter into the sheath advancing it until all the fenestrations are within the pleural cavity. This can be verified under fluoroscopy as fenestrations are located along the barium sulfate stripe.
18. Peel away the sheath while ensuring the catheter remains in place. Adjust the catheter so that it lies flat in the tunnel without any kinks.

**Caution:** Do not use forceps on the introducer to break the handle and/or peel the sheath.

19. Close the incision at the guidewire insertion site.
20. Close the incision site around the catheter and suture the catheter to the skin taking care not to restrict the diameter of the catheter. This suture is intended to remain in place at least until there is tissue ingrowth around the cuff.

**Caution:** Exercise care when placing ligatures to avoid cutting or occluding the catheter.

**Note:** After suturing, Dermabond™ Topical Skin Adhesive (not included) may be applied topically over the guidewire insertion site and catheter exit site. This may aid in preventing leakage by assisting in wound closure. Follow the instructions for use provided by the manufacturer.

#### Drainage Procedure

The drainage procedure can be performed using:

- a) PleurX Vacuum Bottle(s)
- b) PleurX Lockable Drainage Line with glass vacuum bottle(s) or with wall suction.

If using PleurX Vacuum Bottle(s), refer to PleurX Drainage Kit Instructions for Use to perform the drainage procedure.

**Caution:** Re-expansion pulmonary edema may occur if too much fluid is removed too rapidly. Therefore, it is recommended to limit the initial drainage to no more than 1,500 ml. The volume of pleural fluid removed should be based on the patient's individual status.

**Caution:** Potential complications of access and drainage of the pleural cavity include, but may not be limited to, the following: re-expansion pulmonary edema, pneumothorax, laceration of the lung or liver, hypotension/circulatory collapse, wound infection, empyema and infection in the pleural cavity.

#### Connect the Drainage Line to Wall Suction

**Caution:** Keep the valve on the PleurX Catheter and the lockable access tip on the drainage line clean. Keep them away from other objects to help avoid contamination.

**Caution:** If wall suction or portable suction is used, it must be regulated to no greater than -60 cm H<sub>2</sub>O, or to drain no more than 400 ml of fluid per minute. (-60 cm H<sub>2</sub>O = -1.7 in Hg = -44 mm Hg = -0.8 psi) Make sure that the portable suction unit is capable of regulating to -60 cm H<sub>2</sub>O or less. Check portable suction unit instructions for use prior to draining.

1. Close the roller clamp completely by rolling the wheel on the roller clamp toward the suction source. **(Figure 6)**

**Caution:** The roller clamp must be fully closed to occlude the drainage line. When not connected to a suction source, make sure the roller clamp is fully closed; otherwise the drainage line may allow air into the body or let fluid leak out.

2. Attach the 5-in-1 adapter to the Luer fitting on the drainage line.
3. Connect the 5-in-1 adapter to the suction source.
4. Remove the cover with the lockable access tip by twisting it and pulling gently. Discard the cover. **(Figure 7)**
5. Insert the lockable access tip on the drainage line securely into the catheter valve. You will feel and hear a click when the lockable access tip and valve are securely connected. **(Figure 8)**
6. If desired, lock the access tip to the catheter valve by twisting the access tip until you feel and hear a second click. **(Figure 9)**

**Caution:** Make sure that the valve and the lockable access tip are securely connected when draining. If they are accidentally separated, they may become contaminated. If this occurs, clean the valve with an alcohol pad and use a new drainage line to avoid potential contamination.

**Caution:** Precautions should be taken to ensure the drainage line is not tugged or pulled.

#### Connect the Drainage Line to Glass Vacuum Bottle(s)

**Caution:** Keep the valve on the PleurX Catheter and the lockable access tip on the drainage line clean. Keep them away from other objects to help avoid contamination.

**Caution:** When draining with glass vacuum bottles, do not use a needle larger than 17 G.

1. Close the roller clamp completely by rolling the wheel on the roller clamp toward the glass vacuum bottle. **(Figure 6)**

**Caution:** The roller clamp must be fully closed to occlude the drainage line. When not connected to a glass vacuum bottle, make sure the roller clamp is fully closed; otherwise the drainage line may allow air into the body or let fluid leak out.

**Caution:** When connecting to a glass vacuum bottle, make sure the roller clamp on the drainage line is fully closed. Otherwise, it is possible for some or all of the vacuum in the bottle to be lost.

2. Attach a 17 G needle to the Luer fitting on the drainage line.
3. Connect the 17 G needle to the glass vacuum bottle.
4. Remove the cover with the lockable access tip by twisting it and pulling gently. Discard the cover. **(Figure 7)**
5. Insert the lockable access tip on the drainage line securely into the catheter valve. You will feel and hear a click when the lockable access tip and valve are securely connected. **(Figure 8)**
6. If desired, lock the access tip to the catheter valve by twisting the access tip until you feel and hear a second click. **(Figure 9)**

**Caution:** Make sure that the valve and the lockable access tip are securely connected when draining. If they are accidentally separated, they may become contaminated. If this occurs, clean the valve with an alcohol pad and use a new drainage line to avoid potential contamination.

**Caution:** Precautions should be taken to ensure the drainage line is not tugged or pulled.

#### Drain Fluid

**Caution:** It is normal for the patient to feel some discomfort or pain when draining fluid. If discomfort or pain is experienced when draining, roll the wheel on the roller clamp towards the suction source to slow or stop the flow of fluid for a few minutes. Pain may be an indication of infection.

1. Roll the wheel on the roller clamp away from the suction source to begin drainage. **(Figure 10)**  
When fluid begins to drain, you may roll the wheel on the roller clamp toward the suction source to slow the flow of fluid.
2. If you need to change the vacuum bottle / canister / water seal device or suction source for any reason, remove the drainage line from the suction source and connect to a new suction source. Roll the wheel on the roller clamp away from the suction source to resume draining.
3. When fluid flow stops or the desired amount of fluid has been removed, completely close the roller clamp by rolling the wheel on the roller clamp toward the suction source. **(Figure 6)**

#### Finish Drainage

1. If locked, twist the lockable access tip to unlock it from the catheter valve. **(Figure 11)**
2. Ensure the drainage line has been unlocked. Pull the lockable access tip out of the valve in a firm, smooth motion. Set the used drainage line down. **(Figure 12)**

3. Clean the catheter valve with an alcohol pad. Do not try to push anything through the valve as damage to the valve may occur. **(Figure 13)**
4. Place the new cap over the catheter valve and twist it clockwise until it clicks into its locked position. Cap might feel loose until it is locked into place. **(Figure 14)**
5. Disconnect the drainage line from the glass vacuum bottle or suction source.

**Note:** Dispose of the used drainage line and/or used vacuum bottles in accordance with applicable local, state, and federal regulations. Used product may present a potential biohazard.

#### Place the Self-Adhesive Dressing

1. Clean around the catheter site.

**Note:** Make sure the exit site and skin around the catheter are dry before completing the dressing procedure.

2. Place the foam catheter pad around the catheter. **(Figure 15)**
3. Wind the catheter into loops and place it over the foam pad. **(Figure 16)**
4. Cover the catheter with up to four (4) gauze pads. **(Figure 17)**
5. The self-adhesive dressing has three (3) layers:
  - a. printed liner
  - b. clear wound dressing
  - c. center panel and frame backing
6. Peel the printed liner from the self-adhesive dressing, exposing the adhesive surface. **(Figure 18)**
7. Center the self-adhesive dressing over the gauze pads and press it down. **(Figure 19)** **Note:** Do not stretch the self-adhesive dressing during application.
8. Remove and discard the center panel from the backing of the self-adhesive dressing. **(Figure 20)**
9. Slowly remove the frame while smoothing down the self-adhesive dressing edges. **(Figure 21)**
10. Smooth the entire self-adhesive dressing from the center toward the edges using firm pressure to enhance adhesion.

**Note:** The Ethicon BioPatch™ (not included) is compatible for use with the PleurX Catheter. If using BioPatch, the disc should be placed on the skin first, before the foam pad and other dressing materials. Follow all instructions for use provided by Ethicon and CareFusion.

**Spontaneous Pleurodesis with the PleurX Catheter**

Patients who drain regularly every day or every other day may achieve pleurodesis. In a multi-center clinical trial, drainage of the effusion at least once every other day resulted in approximately half of the patients achieving spontaneous pleurodesis with a mean time to catheter removal of 29 days. <sup>1</sup>

**Pleurodesis Procedures with the PleurX Catheter**

The PleurX Catheter is compatible with Talc and Bleomycin (not included).

**Talc Slurry / Bleomycin**

Use the PleurX Catheter Access Kit (REF 50-7280) to access the catheter and instill the pleurodesing agent. The locking access tip with needleless access valve should be removed immediately after completion of the procedure.

Perform the procedure per the pleurodesis agent's instructions for use and your institution protocol.

**Talc Poudrage**

The PleurX catheter can be used in place of the chest tube for post-surgical drainage of a talc poudrage.

**Note:** The poudrage should not be sprayed through the PleurX catheter. Perform the procedure per the Talc instructions for use and your institution protocol.

**Subsequent Drainage Procedures**

Subsequent drainage procedures are to be performed using the PleurX Lockable Drainage Line, PleurX Vacuum Bottle, or the PleurX Drainage Kits. Each drainage kit contains the necessary drainage line, vacuum bottle, and other necessary items to perform the drainage procedure.

**It is vital that patients and/or caregivers are carefully instructed on how to use the kit to drain the pleural cavity. The person(s) responsible for drainage must be able to demonstrate they are capable of performing the procedure.**

**If the patient/caregiver is not able or willing to perform the drainage, a medical professional should perform the drainage.**

**It is recommended that the patient is periodically contacted or seen by a clinician to evaluate treatment regimen and evaluate catheter functional status.**

**Catheter Maintenance**

To perform saline flushing or dec clotting procedure, use the PleurX Catheter Access Kit (REF 50-7280). If the catheter valve or tubing is damaged, use the PleurX Valve Kit (REF 50-7270) to cut the tubing and replace the valve.

**Catheter Removal Procedure**

It may be appropriate and/or necessary at a later date to remove the PleurX Pleural Catheter. Three successive attempts to drain fluid that result in less than 50 ml of fluid removed may indicate one of the following:

- pleurodesis has been achieved
  - the catheter is loculated away from the fluid
  - the catheter is occluded
1. Place the patient appropriately to access the catheter insertion site.
  2. Aseptically clean the patient's chest around the catheter insertion site.
  3. Anesthetize the site.
  4. Remove any remaining sutures securing the catheter.
  5. Using forceps, dissect around the cuff to free it from the ingrowth. Ensure that the cuff is completely free within the tunnel.
  6. Grasp the catheter in one hand and pull with a firm, constant pressure.
  7. Cover the site as appropriate.

Contains Phthalates. The benefit of treatment outweighs the remote possibility of exposure to phthalates.

**Note:** Not made with natural rubber latex.

**Magnetic Resonance Imaging (MRI) Safety Information**

Non-clinical testing has demonstrated that the PleurX Catheter is MR safe. In non-clinical testing, the image artifact caused by the device extends approximately 2 mm from the device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

<sup>1</sup>Putnam JB Jr, Light RW, Rodriguez RM, et al. A Randomized Comparison of Indwelling Pleural Catheter and Doxycycline Pleurodesis in the Management of Malignant Pleural Effusions. Cancer 1999; 86; 1992-1999.

Proofed by: \_\_\_\_\_ Date: \_\_\_\_\_  
Dimensions checked: \_\_\_\_\_ Copy checked: \_\_\_\_\_

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Becton Dickinson İth. İhr. Ltd.