**REF** PIG1260K 6Fr/PIG1280K 8Fr

## Safe-T-Centesis™

- en** **Safe-T-Centesis Catheter Drainage Kit**  
Do not use if package is damaged.
- fr** **Kit de drainage pour cathéter Safe-T-Centesis**  
Ne pas utiliser si l'emballage est endommagé.
- de** **Safe-T-Centesis-Katheterdrainagekit**  
Nicht verwenden, wenn die Verpackung beschädigt ist.
- it** **Kit di drenaggio per catetere Safe-T-Centesis**  
Non usare se la confezione è danneggiata.
- es** **Juego de drenaje de catéter Safe-T-Centesis**  
No utilizar si el paquete está dañado.
- pt** **Kit de cateter de drenagem Safe-T-Centesis**  
Não utilizar se a embalagem estiver danificada.
- nl** **Drainagekit Safe-T-Centesis-katheters**  
Niet gebruiken indien de verpakking beschadigd is.
- sv** **Safe-T-Centesis-dräneringskit för kateter**  
Får ej användas om förpackningen är skadad.
- da** **Dræningssæt med Safe-T-Centesis-kateter**  
Anvend ikke, hvis emballagen er beskadiget.
- fi** **Katetrin Safe-T-Centesis -dreneerausarja**  
Ei saa käyttää, jos pakkaus on vaurioitunut.
- no** **Safe-T-Centesis kateterdreneringssett**  
Produktet må ikke brukes hvis forpakningen er skadd.
- el** **Κιτ παροχέτευσης καθετήρα Safe-T-Centesis**  
Μην το χρησιμοποιείτε, εάν η συσκευασία έχει φθαρεί.
- ja** **Safe-T-Centesis カテーテルドレナージキット**  
包装が開封されるか又は損傷を受けている場合は使用しないでください。

Not made  
with natural  
rubber latex.

**USA** Rx Only

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
## Instructions for Use

## Caution

- Ensure all Luer-Lock connections are secure before using device.
- While the blunt obturator and color change indicator have been incorporated to minimize the likelihood of a penetrating injury, standard precautionary measures employed in all needle insertions must be observed.
- Do not attempt patient insertion if the color change indicator does not transition (white/red) after complete mating of the safety introducer needle assembly and pigtail catheter assembly.
- Ensure subcutaneous placement of all catheter ports prior to safety introducer needle removal.
- Upon completion of the drainage procedure, remove the catheter quickly to avoid entry of air.
- This device is not intended for use with sclerosing or chemotherapeutic agents.
- This device is intended for single patient use.

**Caution:** For single use only. Re-use may result in a non-functional product or contribute to cross contamination.

Do not use if packaged is damaged. Do not resterilize.

- A. Safety Introducer Needle Assembly
- B. Pigtail Catheter Assembly
- C. Color Change Indicator
- D. Threaded Valve Assembly
- E. Remote Stopcock
- F. Safe-T-Centesis Device (A+B=F)

## Procedure

After puncture site has been prepared and anesthetized:

1. The safety scalpel can be used to make a small incision.
2. Remove both pieces of the needle protector from the safety introducer needle assembly.
3. Gently load the safety introducer needle assembly into the proximal end of the pigtail catheter drainage assembly. Ensure a complete engagement of the threaded valve assembly and the safety introducer needle housing assembly to unlock the blunt obturator.
4. Ensure the sharp portion of the safety introducer needle assembly extends beyond the catheter tip.
5. Introduce the Safe-T-Centesis Device for either: **Thoracentesis** into the chest wall or **Paracentesis** into the abdominal region utilizing sufficient pressure to detect movement of the color change indicator within the safety introducer needle housing assembly. The indicator should change from white to red when the blunt obturator is retracted into the cannula and the sharp needle is exposed.
6. Advance the Safe-T-Centesis Device into the desired cavity. The color change indicator should transition from red to white and an audible sound may be heard as the needle enters a "free cavity." Failure of the color change indicator to transition back to white indicates the sharp needle tip is exposed. If the function of the device is in question, the device should be discarded.
7. Utilize the attached syringe to aspirate fluid and confirm positive placement.
8. Ensure the stopcock flow control arm is turned off to pigtail drainage assembly and advance the pigtail drainage assembly off of the safety introducer needle assembly into effusion in accordance with physician's clinical judgment. **Caution: Reinsertion of the safety introducer needle assembly after full removal from the pigtail catheter drainage assembly may compromise the integrity of the valve seal.**
9. **These steps apply when a syringe is used for fluid collection:**
  - A. Attach a syringe to back port of the stopcock using a push-twist motion.
  - B. Turn the stopcock off to the side port.
  - C. Collect fluid.

**These steps apply when a vacuum bottle is used for fluid collection:**

  - A. Attach male Luer end of extension set to rear port of stopcock.
  - B. Tighten all connection with push-twist motion prior to proceeding.
  - C. Remove protector from plastic vacuum bottle spike and puncture vacuum bottle. If desired, the additional 16 G x 1 1/2" vacuum bottle needle can be used to puncture vacuum bottle: remove plastic vacuum bottle spike from extension set, replace with 16 G x 1 1/2" vacuum bottle needle, remove protector and puncture vacuum bottle.
  - D. Turn stopcock off to side port. This will start fluid flow into the bottle.

**These steps apply when wall suction is used:**

  - A. Attach male Luer end of extension set to rear port of stopcock.
  - B. Remove and dispose of plastic vacuum bottle spike with protector and replace with 5-in-1 adapter.
  - C. Tighten all connections with a push-twist motion prior to proceeding.
  - D. Insert 5-in-1 adapter firmly into wall suction drainage tubing.
  - E. Start wall suction drainage and turn stopcock off to side port.
10. Upon completion of procedure, turn the stopcock off to the pigtail drainage assembly.
11. Remove the drainage catheter from the patient and protect the puncture site with a bandage. If the additional 16 G x 1 1/2" vacuum bottle needle has been utilized, the Point-Lok needle retention device should be utilized as a sharps engineered safety device.

These instructions are supplied to explain the use of components in this kit. Procedural details will vary according to the physician's clinical judgment. Contains Phthalates. The benefit of treatment outweighs the remote possibility of exposure to phthalates.

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0123

Made in Dominican Republic by  
CareFusion DR 203 LTD  
362-23401 • 2017-11

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© 2010, 2012, 2014, 2017 CareFusion Corporation or one of its subsidiaries. All rights reserved.

Safe-T-Centesis, CareFusion and the CareFusion logo are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries.

Point-Lok is a registered trademark of Smiths Medical ASD.

Sponsored by  
CareFusion Australia 316 Pty. Ltd.  
CareFusion New Zealand 313 Ltd.

DEHP