REF 50-9050A

PeritX™ Peritoneal Catheter Mini Kit

- en** Peritoneal Catheter Mini Kit
For Peritoneal Placement Only
- fr** Mini kit de cathéter péritonéal
Pour un placement péritonéal uniquement
- de** Peritonealkatheter-Minikit
Nur zur peritonealen Anwendung
- it** Mini kit per catetere peritoneale
Esclusivamente per posizionamento peritoneale
- es** Mini kit de catéter peritoneal
Solo para colocación peritoneal
- pt** Mini kit de cateter peritoneal
Apenas para colocação peritoneal
- nl** Minikit voor peritoneale katheter
Uitsluitend voor peritoneale plaatsing
- sv** Peritonealkateter, miniset
Endast för peritoneal placering
- fi** Peritoneaalikatetriminiseti
vain vatsakalvoon asetettavaksi
- no** Peritonealkateter – miniset
Bare for peritoneal plassering
- da** Peritonealkateter, minisæt
Kun til anlæggelse i peritoneum
- el** Μίνι κιτ περιτοναϊκού καθετήρα
Για περιτοναϊκή τοποθέτηση μόνο
- tr** Periton Kateteri Mini Kit
Yalnız Periton Boşluğuna Yerleştirilmeye Yöneliktir

Not made with
natural rubber
latex.

STERILE **EO**

For Single Use Only

Do Not Resterilize

Do not use if
package is damaged.

DEHP

Contains phthalates

CareFusion
75 North Fairway Drive
Vernon Hills, IL 60061 USA

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Becton Dickinson Pty Ltd
Level 5, 66 Waterloo Road
Macquarie Park, NSW 2113, Australia

Becton Dickinson Limited
14B George Bourke Drive
Mt. Wellington, Auckland, 1060 New Zealand
361-42401 • 2020-11

USA Rx Only

© 2020 CareFusion Corporation or one of its affiliates. All rights reserved.

For product inquiries: bd.com/drainage

PleurX, PeritX, CareFusion and the CareFusion logo are trademarks or registered trademarks of CareFusion Corporation, or one of its affiliates.

BioPatch and Dermabond are trademarks of Johnson & Johnson Corporation.


CareFusion


Figure, Abbildung, Figura, Afbeelding, Figur, Kuva, Εικόνα, Şekil


Figure, Abbildung, Figura, Afbeelding, Figur, Kuva, Eικόνα, Şekil

(4)


(5)

(6)


"CLICK" "CLIQUEZ" "KLICK"
"CLIC" "KLIK" "CLIQUE"

"KLIKK PÅ" "NAKSAHDUS"
"KLIK" "κλικ" "klik"


(8)


"CLICK" "CLIQUEZ" "KLICK"
"CLIC" "KLIK" "CLIQUE"

"KLIKK PÅ" "NAKSAHDUS"
"KLIK" "κλικ" "klik"


(9)

(10)









(7)

(11)

Figure, Abbildung, Figura, Afbeelding, Figur, Kuva, Εικόνα, Şekil


PeritX™ Peritoneal Catheter Procedure


Product Information

The PeritX™ Peritoneal Catheter consists of a fenestrated silicone catheter with a valve mechanism and a polyester cuff. A barium sulfate stripe runs the entire length of the catheter. (Figure 1) The valve is designed to prevent the passage of air or fluid in either direction unless it is accessed with the specifically matched drainage line, access kit, or vacuum bottles provided by CareFusion. The PeritX™ Peritoneal Catheter is designed exclusively for use with the PleurX™ Vacuum Bottles, glass vacuum bottles, and the Lockable Drainage Line for connection to wall suction or portable suction.

PeritX™ Peritoneal Catheter

Lockable Drainage Line

Indications for Use

The PeritX™ Peritoneal Catheter Mini Kit System is indicated for intermittent, long term drainage of symptomatic, recurrent, malignant and non-malignant ascites that does not respond to medical management of the underlying disease and for the palliation of symptoms related to recurrent ascites. The use of the PeritX™ Peritoneal Catheter for non-malignant ascites is limited to patients who are intolerant or resistant to maximum medical therapy, refractory to large volume paracentesis (LVP) and are not candidates for a trans-jugular intrahepatic portosystemic shunt or LVP. The PeritX™ Peritoneal Catheter is indicated for adults only.

The Lockable Drainage Line (Figure 2) is used to drain fluid using standard wall suction, water seal drainage system, vacuum bottle or other appropriate method.

Contraindications

Use of the PeritX™ Peritoneal Catheter System is contraindicated in the following situations:

1. When the peritoneal cavity is multi-loculated, and the drainage of a single loculation would not be expected to provide relief of dyspnea or other symptoms.

2. When there is a coagulopathy.
3. When the peritoneal cavity is infected.
4. When the patient is known or suspected to be allergic to materials contained in the device.

Warnings

Do not put anything except the access tip of the Lockable Drainage Line, Catheter Access Kit, or PleurX™ Vacuum Bottles into the catheter valve since any other device could damage the valve. A damaged valve may allow air into the body or let fluid leak out through the valve when not draining.

A diagnostic paracentesis should be performed if the patient shows signs or symptoms of possible spontaneous bacterial peritonitis (SBP) such as fever or abdominal pain. If SBP is present, the patient should be treated per institutional guidelines, including systemic antibiotics and repeat diagnostic paracentesis at the end of the antibiotic regimen. For a patient with resolved SBP, the patient should be treated per institutional guidelines, including prophylactic antibiotics to help in the prevention of refractory or recurrent SBP. In the case of refractory or recurrent infections, the catheter should be removed and reinserted at the discretion of the clinician after the SBP has resolved.

Do not place the PeritX™ Peritoneal Catheter into the pleural space as it could lead to misidentification of the catheter and/or mistreatment of the patient.

Cautions

For single use only. Re-use may result in a non-functional product or contribute to cross contamination.

Sterile technique should be used when placing and draining the catheter.

Care must be taken when inserting the guidewire introducer needle to avoid puncturing or lacerating the liver or bowel.

If the needle is left in place with the guidewire introducer, damage to the guidewire may result if withdrawn through the needle.

The fenestrations must be entirely within the peritoneal space to avoid leakage into the tunnel tract. Take patient size, tunnel length and the catheter length into account when placing the catheter.

Potential complications of access and drainage of the peritoneal cavity include, but may not be limited to, the following: laceration of liver or bowel, hypotension/circulatory collapse, electrolyte imbalance, protein depletion, ascites leakage, peritonitis, wound infection, intraperitoneal adhesion, tumor growth in the catheter tunnel, and loculation of the peritoneal cavity.

Removal of chylous ascites could exacerbate protein depletion or related nutritional complications.

Sterility

This product has been sterilized. It is for single use only and is not to be resterilized. Do not use if package is damaged. CareFusion will not be responsible for any product that is resterilized, nor accept for credit or exchange, any product that has been opened but not used.

PeritX™ Catheter Tray Supplies

Placement Components

- 1 PeritX™ Peritoneal Catheter, 15.5 Fr
- 1 Guidewire Introducer with Needle
- 1 Syringe, 10 ml
- 1 J-Tip Guidewire
- 1 Peel-Away Introducer, 16 Fr
- 1 Tunneler
- 1 Needle Foam Stop
- 1 Blue Wrap

Drainage Components

- 1 Drainage Line with Lockable Access Tip
- 1 Needle, 17 G x 1" (2.5 cm)
- 1 5-in-1 Drainage Line Adapter
- 1 Valve Cap

Dressing Components


- 6 Gauze Pads, 4" x 4" (10.2 cm x 10.2 cm)
- 1 Foam Catheter Pad
- 1 Self-Adhesive Dressing

General Guidelines

1. Systemic prophylactic antibiotics may be indicated.
2. The procedure for peritoneal placement can be performed using local anesthetic and sedation. However, depending on patient needs, it may be performed using alternative approaches to anesthesia or sedation.
3. The catheter should be placed under image guidance, using all precautions normally used for percutaneous placement of indwelling, tunneled catheters. Care should be taken to identify and avoid contact with vasculature near the guidewire insertion site.
4. Guidewire insertion site selection should be based upon patient anatomy and presentation with consideration given to any possible adhesions or loculated pockets of fluid. The fenestrated section of the catheter should preferentially be placed low in the peritoneal cavity to maximize access to fluid.
5. Consideration should be given to the patient's ease of access in determining the location of the catheter exit site.
6. The fenestrated end of the catheter may be cut shorter depending on an individual patient's anatomy. **Note:** Leave at least one fenestration present on the catheter.
7. When using the drainage line, it may be used to drain at a suction level of -60 cm H₂O or less for as long as it takes to drain 2,000 ml of fluid from the abdomen per day. The volume of fluid removed should be based on the individual patient's status and the risks for over draining, including hypotension.

8. The use of this product for refractory non-malignant ascites should be limited to patients who have failed 1 month of maximum medical therapy and 3 LVP procedures within 3 months.
9. In order to ensure the catheter is not placed in a patient with a pre-existing infection, a diagnostic paracentesis evaluating for evidence of infection (such as gram stain and/or Absolute Neutrophil Count) should be performed prior to the placement of the device.
10. Drainage of the peritoneal cavity prior to catheter placement should be performed in the case of tense ascites to provide acute relief and reduce abdominal wall tension. Complete drainage prior to catheter placement is not recommended.

Suggested Placement Procedure

Proper medical and surgical procedures are the responsibility of the physician. The appropriateness of any procedure must be based upon the needs of the patient. **Figure (3)** illustrates the placement of the PeritX™ Peritoneal Catheter, as described in the following procedure.

Caution: Individual patient anatomy, such as thin or weak abdominal wall, may require procedural variations to reduce the risk of leakage around the catheter.

1. Position the patient appropriately to access the desired guidewire insertion site.
2. Identify the appropriate insertion site through which to place the guidewire. The guidewire insertion site is typically lateral to the midline, 6–10 cm below the costal margin, and above the patient's beltline. Ultrasound can be used to confirm the guidewire insertion site.
3. Identify the location of the catheter exit site, which is usually 5-8 cm from the guidewire insertion site.
4. Surgically prep the patient.
5. Drape and anesthetize the planned insertion and tunneling sites.

Note: It is important to use an oblique approach while passing the introducer needle through the peritoneal wall to minimize the chance of ascitic fluid leakage.

Caution: Care must be taken when inserting the needle to avoid puncturing or lacerating the liver or bowel.

6. Insert the guidewire introducer with needle, attached to a syringe, through the peritoneum using an oblique angle.
7. Ensure free aspiration of ascitic fluid, then; remove the needle and syringe, leaving the guidewire introducer in place.
8. Insert the guidewire through the introducer and advance the distal end into the peritoneal cavity.

Caution: Do not allow the guidewire to inadvertently advance totally within the patient. Ensure that the guidewire exits the proximal end of the needle, dilator, or sheath prior to and during placement of each component.

9. Remove the introducer, leaving the guidewire in place.

Caution: Damage to the guidewire may result if withdrawn through the needle.

10. Make a 1 cm incision at the guidewire insertion site.
11. Make a second 1-2 cm incision 5-8 cm from the guidewire insertion site. This incision will be the catheter exit site. Consider the patient's ease of access in determining its location. **(Figure 3)**

Note: A smaller incision may provide better security of the catheter.

Note: Take care to ensure that the tunnel track has been anesthetized.

12. The fenestrated end of the catheter may be cut shorter depending on an individual patient's anatomy. If desired, cut a portion of the fenestrated end of the catheter by using a scalpel to make a straight cut between fenestrations. Leave at least one (1) fenestration on the catheter.
13. Attach the fenestrated end of the catheter onto the tunneler.

Caution: Exercise care when handling the catheter to prevent it from coming into contact with surfaces such as drapes or towels. Silicone rubber is highly electrostatic and attracts airborne particles and surface contaminants.

Caution: Use smooth surface instruments when handling the catheter. Possible cuts or tears can occur if smooth surface instruments are not used.

14. Pass the tunneler **(A)** and catheter **(B)** subcutaneously from the second incision down to and out through the first incision at the guidewire insertion site. **(Figure 4)** Continue to draw the catheter through the tunnel until the polyester cuff lies inside the tunnel, about 1 cm **(C)** from the second incision. **(Figure 5)** Disconnect the tunneler from the catheter.

Note: If the cuff is advanced further into the tunnel, it can make later removal of the catheter difficult.

15. Thread the 16 Fr peel-away introducer sheath over the guidewire into the peritoneal cavity.
16. Remove the guidewire and dilator as a unit, leaving the 16 Fr peel-away introducer sheath in place.

Caution: Place a thumb over the end of the sheath as the dilator is removed to avoid peritoneal fluid leakage. Care must be taken not to bend or kink the sheath. Damage to the sheath may prevent passage of the catheter.

17. Insert the fenestrated end of the catheter into the sheath advancing it until all the fenestrations are within the peritoneal cavity. This can be verified under fluoroscopy as fenestrations are located along the barium sulfate stripe.
18. Peel away the sheath while ensuring the catheter remains in place. Adjust the catheter so that it lies flat in the tunnel without any kinks.

Caution: Do not use forceps on the introducer to break its handle and/or peel the sheath.

19. Close the incision at the guidewire insertion site.
20. Close the incision site around the catheter and suture the catheter to the skin taking care not to restrict the diameter of the catheter. This suture is intended to remain in place at least until there is tissue ingrowth around the cuff.

Caution: Exercise care when placing ligatures to avoid cutting or occluding the catheter.

Note: After suturing, Dermabond™ Topical Skin Adhesive (not included) may be applied topically over the guidewire insertion site and catheter exit site. This may aid in preventing leakage by assisting in wound closure. Follow the instructions for use provided by the manufacturer.

Drainage Procedure

The drainage procedure can be performed using:

- a) PleurX™ Vacuum Bottle(s)
- b) Lockable Drainage Line with glass vacuum bottle(s) or with Wall Suction

If using PleurX™ Vacuum Bottle(s), refer to PleurX™ Drainage Kit Instructions for Use.

Caution: In ascites patients, paracentesis-related hypotension is uncommon, but has been documented. Use of IV fluid replacement and/or administration of colloidal agents can reduce the risk of hypotension. Additionally, initial drainage should be no more than 6L in the first 24 hours.

Connect the Drainage Line to Wall Suction

Caution: Keep the valve on the catheter and the lockable access tip on the drainage line clean. Keep them away from other objects to help avoid contamination.

Caution: If wall suction is used, it must be regulated to no greater than -60 cm H₂O, or to drain no more than 400 ml of fluid per minute.
(-60 cm H₂O = -1.7 in Hg = -44 mm Hg = -0.8 psi)

1. Close the roller clamp completely by rolling the wheel on the roller clamp toward the suction source. **(Figure 6)**

Caution: The roller clamp must be fully closed to occlude the drainage line. When not connected to a suction source, make sure the roller clamp is fully closed; otherwise the drainage line may allow air into the body or let fluid leak out.

2. Attach the 5-in-1 adapter to the Luer fitting on the drainage line.
3. Connect the 5-in-1 adapter to the suction source.
4. Remove the cover with the lockable access tip by twisting it and pulling gently. Discard the cover. **(Figure 7)**
5. Insert the lockable access tip on the drainage line securely into the catheter valve. You will feel and hear a click when the lockable access tip and valve are securely connected. **(Figure 8)**
6. If desired, lock the access tip to the catheter valve by twisting the access tip until you feel and hear a second click. **(Figure 9)**

Caution: Make sure that the valve and the lockable access tip are securely connected when draining. If they are accidentally separated, they may become contaminated. If this occurs, clean the valve with an alcohol pad and use a new drainage line to avoid potential contamination.

Caution: Precautions should be taken to ensure the drainage line is not tugged or pulled.

Connect the Drainage Line to Glass Vacuum Bottle(s)

Caution: Keep the valve on the PeritX™ Catheter and the lockable access tip on the drainage line clean. Keep them away from other objects to help avoid contamination.

Caution: When draining with glass vacuum bottles, do not use a needle larger than 17 G.

1. Close the roller clamp completely by rolling the wheel on the roller clamp toward the glass vacuum bottle. **(Figure 6)**

Caution: The roller clamp must be fully closed to occlude the drainage line. When not connected to a glass vacuum bottle, make sure the roller clamp is fully closed; otherwise the drainage line may allow air into the body or let fluid leak out.

Caution: Before connecting to a glass vacuum bottle, close the roller clamp completely by rolling the wheel on the roller clamp towards the glass vacuum bottle; otherwise it is possible for some, or all of the vacuum in the bottle to be lost.

2. Attach a 17 G needle to the Luer fitting on the drainage line.
3. Connect the 17 G needle to the glass vacuum bottle.
4. Remove the cover with the lockable access tip by twisting it and pulling gently. Discard the cover. **(Figure 7)**

5. Insert the lockable access tip on the drainage line securely into the catheter valve. You will feel and hear a click when the lockable access tip and valve are securely connected. **(Figure 8)**
6. If desired, lock the access tip to the catheter valve by twisting the access tip until you feel and hear a second click. **(Figure 9)**

Caution: Make sure that the valve and the lockable access tip are securely connected when draining. If they are accidentally separated, they may become contaminated. If this occurs, clean the valve with an alcohol pad and use a new drainage line to avoid potential contamination.

Caution: Precautions should be taken to ensure the drainage line is not tugged or pulled.

Drain Fluid

Caution: It is normal for the patient to feel some discomfort or pain when draining fluid. If discomfort or pain is experienced when draining, roll the wheel on the roller clamp towards the suction source to slow or stop the flow of fluid for a few minutes. Pain may be an indication of infection.

1. Roll the wheel on the roller clamp away from the suction source to begin drainage. **(Figure 10)**
When fluid begins to drain, you may roll the wheel on the roller clamp toward the suction source to slow the flow of fluid.
2. If you need to change the glass vacuum bottle/canister/water seal device or suction source for any reason, remove the drainage line from the suction source and connect to a new suction source. Roll the wheel on the roller clamp away from the suction source to resume draining.
3. When fluid flow stops or the desired amount of fluid has been removed, completely close the roller clamp by rolling the wheel on the roller clamp toward the suction source. **(Figure 6)**

Finish Drainage

1. If locked, twist the lockable access tip to unlock it from the catheter valve. **(Figure 11)**
2. Ensure the drainage line has been unlocked. Pull the lockable access tip out of the valve in a firm, smooth motion. Set the used drainage line down. **(Figure 12)**
3. Clean the catheter valve with an alcohol pad. Do not try to push anything through the valve as damage to the valve may occur. **(Figure 13)**
4. Place the new cap over the catheter valve and twist it clockwise until it clicks into its locked position. Cap might feel loose until it is locked into place. **(Figure 14)**
5. Disconnect the drainage line from the glass vacuum bottle or suction source.

Note: Dispose of the used drainage line and/or used glass vacuum bottles in accordance with applicable local, state, and federal regulations. Used product may present a potential biohazard.

Place the Self-Adhesive Dressing

1. Clean around the catheter site.

Note: Make sure the exit site and skin around the catheter are dry before completing the dressing procedure.

2. Place the foam catheter pad around the catheter. **(Figure 15)**
3. Wind the catheter into loops and place it over the foam pad. **(Figure 16)**
4. Cover the catheter with up to four (4) gauze pads **(Figure 17)**.
5. The self-adhesive dressing has three (3) layers:
 - a. printed liner
 - b. clear wound dressing
 - c. center panel and frame backing
6. Peel the printed liner from the self-adhesive dressing, exposing the adhesive surface. **(Figure 18)**
7. Center the self-adhesive dressing over the gauze pads and press it down. **(Figure 19)** **Note:** Do not stretch the self-adhesive dressing during application.
8. Remove and discard the center panel from the backing of the self-adhesive dressing. **(Figure 20)**
9. Slowly remove the frame while smoothing down the self-adhesive dressing edges. **(Figure 21)**
10. Smooth the entire self-adhesive dressing from the center toward the edges using firm pressure to enhance adhesion.

Note: The Ethicon BioPatch™ (not included) is compatible for use with the PeritX™ Catheter. If using BioPatch, the disc should be placed on the skin first, before the foam pad and other dressing materials. Follow all instructions for use provided by Ethicon and CareFusion.

It is vital that patients and/or caregivers are carefully instructed on how to use the kit to drain ascites. The person(s) responsible for drainage must be able to demonstrate they are capable of performing the procedure.

If the patient/caregiver is not able or willing to perform the drainage, a medical professional should perform the drainage.

It is recommended that the patient is periodically contacted or seen by a clinician to evaluate treatment regimen, assess need for albumin supplementation and evaluate catheter's functional status.

Subsequent Drainage Procedures

Subsequent drainage procedures are to be performed using the Lockable Drainage Line, PleurX™ Vacuum Bottle, or the PleurX™ Drainage Kits. Each drainage kit contains the necessary drainage line, vacuum bottle, and other necessary items to perform the drainage procedure.

Catheter Maintenance

To perform saline flushing or declotting procedure, use the Catheter Access Kit (REF 50-7280A). If the catheter valve or tubing is damaged, use the PeritX™ Valve Kit (REF 50-7272) to cut the tubing and replace the valve.

Catheter Removal Procedure

It may be appropriate and/or necessary at a later date to remove the PeritX™ Peritoneal Catheter. Three successive attempts to drain fluid that result in less than 50 ml of fluid removed may indicate one of the following:

- the ascites has resolved
 - the catheter is loculated away from the fluid
 - the catheter is occluded
1. Place the patient appropriately to access the catheter insertion site.
 2. Aseptically clean the patient's abdomen around the catheter insertion site.
 3. Anesthetize the site.
 4. Remove any remaining sutures securing the catheter.
 5. Using forceps, dissect around the cuff to free it from the ingrowth. Ensure that the cuff is completely free within the tunnel.
 6. Grasp the catheter in one hand and pull with a firm, constant pressure.
 7. Cover the site as appropriate.

Caution: Additional complications during catheter removal procedure may be possible due to long catheter indwell time.

Contains Phthalates. The benefit of treatment outweighs the remote possibility of exposure to phthalates.

Note: Not made with natural rubber latex.

Magnetic Resonance Imaging (MRI) Safety Information

Non-clinical testing has demonstrated that the PeritX™ Catheter is MR safe. In non-clinical testing, the image artifact caused by the device extends approximately 2 mm from the device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

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Becton Dickinson İth. İhr. Ltd.