ewimed

DE drainova® Katheter-Set / drainova® ArgentiC Katheter-Set

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EN drainova® catheter kit / drainova® ArgentiC catheter kit

Instructions for use

FR drainova® kit cathéter / drainova® ArgentiC kit cathéter

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EN drainova® catheter kit/ drainova® ArgentiC catheter kit

- | | |
|--|---|
| 1 drainova® catheter/ drainova® ArgentiC catheter | 11 Foam pad |
| 2 drainova® tunneler | 12 Self-adhesive film dressing |
| 3 Tearaway introducer 16 Fr | 13 Emergency clamp |
| 4 Introducer needle 18 G | 14 drainova® silicone cap, sterile |
| 5 Syringe 10 ml | 15 Surgical drape |
| 6 Guidewire with J-Tip | 16 Surgical drape with aperture |
| 7 drainova® flushing line | 17 Surgical scissors |
| 8 drainova® cone adapter | 18 Mayo-Hegar needle holder |
| 9 drainova® clickFix | 19 Chiraflon suture |
| 10 Gauze pads | 20 Scalpel no. 11P |

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CONTENT

2	Explanation of symbols	36	Recommended Implantation Procedure
7	Product description		Pleural Effusion
30	drainova® Catheter Kit / drainova® ArgentiC Catheter Kit	38	Recommended Implantation Procedure
30	Intended Purpose		Ascites
31	The drainova® Catheter	41	Drainage Procedure
31	The drainova® ArgentiC Catheter	42	Completion of Drainage Procedure
31	Warnings	42	Connecting the Rinsing Tube and Flushing the Catheter (Dislodging Occlusions)
	general	43	Catheter Explanation
33	Warnings	45	Glossary
	regarding implantation		
34	Warnings		
	post implantation		
35	General Preparation		

EN

Please carefully read the instructions for use before using the drainova® Catheter Kit / drainova® ArgentiC Catheter Kit.

The drainova® Catheter Kit and drainova® ArgentiC Catheter Kit will be hereinafter referred to as ‚catheter kit‘.

All explanations of terms (*) can be found in the glossary on page 45.

YOU HAVE QUESTIONS OR NEED HELP?

Contact us:

DE +49 7471 73972-0
info@ewimed.com

SE +46 8 25 11 69
info.se@ewimed.com

BE +49 7471 73972-0
info@ewimed.com

EN +49 7471 73972-0
info@ewimed.com

DA +45 31 43 50 60
info.dk@ewimed.com

NL +49 7471 73972-0
info@ewimed.com

CH +41 31 958 66 66
info.ch@ewimed.com

NO +46 8 25 11 69
info.no@ewimed.com

LU +49 7471 73972-0
info@ewimed.com

AT +43 2231 22500
info@ewimed.com

HU +36 1 998 8295
info@ewimed.hu

For further detailed information please visit our website www.ewimed.com

DRAINOVA® CATHETER KIT / DRAINOVA® ARGENTIC CATHETER KIT

The drainova® Catheter Kit and drainova® ArgentiC Catheter Kit contain all components required for the implantation of a drainova® Catheter / drainova® ArgentiC Catheter used in cases of recurring, therapy-resistant, malignant and non-malignant effusions in serous body cavities. Drainage of the accumulated fluid provides relief from symptoms arising from an effusion.

Medical Indication

The use of the catheter is indicated in the following circumstances:

- Malignant pleural effusion*
- Non-malignant pleural effusion (e.g. heart failure)
- Malignant ascites*
- Non-malignant ascites (e.g. liver failure)

Contraindication

The use of the catheter is contraindicated in the following circumstances:

- The body cavity has several chambers (septa), indicating that no alleviation of dyspnoea, or other symptoms is to be expected after drainage.

- Coagulopathy
- Infection in the body cavity (e.g. pleuritis/peritonitis)
- Lymphatic effusion
- Shift of the mediastinum by more than 2cm to the ipsilateral side of the pleural effusion
- In case of allergy to the materials used in the product

The use of catheter kit components is contraindicated in the following situations:

- Syringe: should not be used other than for the intended purpose. The syringe should not be used for high viscous fluid.
- Guide wire: Since guide wires are used as accessories for the applications mentioned above, the contraindications for these applications also apply:
 - Applications to the central nervous system,
 - Applications to the central circulatory system.
- Tearaway sheath introducer: There is increased risk of pneumothorax for the patient who has severe chronic lung disease. Poor healing may result in the patient who has had irradiation to the anterior chest.

INTENDED PURPOSE

The drainova® catheter kit/ drainova® ArgentiC catheter kit serves as a comprehensive compilation of the components required for the implantation of a drainova® catheter / drainova® ArgentiC catheter used in the treatment of recurrent, therapy refractory, malignant and non-malignant effusion accumulations in serous body cavities to alleviate the symptoms by draining the fluids.

Including for example:

- Malignant pleural effusion
- Non-malignant pleural effusion
- Malignant ascites
- Non-malignant ascites

All set components are CE labelled and sterilized with ethylene oxide gas

THE DRAINOVA® CATHETER

The catheter is made of biocompatible silicone and is 66 cm long. The diameter is 15.5 Fr. The catheter has opposite holes over a distance of 26 cm from one end. This part is implanted and the openings enable the flow of fluid into the catheter. The catheter has a safety valve at the non-implanted end. A polyester sleeve is fitted between the perforated side and the safety valve. The polyester cuff secures the catheter in the subcutaneous* tunnel with the skin tissue. It also


reduces the penetration of germs afterwards. The safety valve closes automatically in the event of a disconnection and prevents the entry of air and the unwanted escape of effusion. The safety valve may only be opened with suitable materials that have the insertion pin specially developed for this purpose. A radiopaque barium sulphate strip runs the entire length of the catheter, marking the position of the catheter in the X-ray image.

EN

THE DRAINOVA® ARGENTIC CATHETER

The shape and handling of the drainova® ArgentiC Catheter is identical to that of the drainova® Catheter described above. It differs only in the material of the tube, which is made out of

biocompatible silicone containing microsilver particles. Silver protects against bacterial colonisation on surfaces and reduces the risk of infection.

WARNINGS

general

Do not use this product if:

- the sterile packaging is damaged
- the product is damaged
- expiry date has elapsed.

The product should solely be used for its intended purpose.

Patient target group: The medical device may only be used in adults.

The drainage consumables are intended for single use only!

Re-use can contribute to a cross-contamination. The catheter is intended for single implantation only.

EN

The catheter kits are delivered in sterile packaging which have undergone quality controls. The catheter kits can not be re-sterilised.

ewimed GmbH accepts no liability for re-used, re-sterilised products or products taken from damaged packaging.

drainova® clickFix: Please make sure to insert the tube connection correctly into the drainova® clickFix to avoid overly squeezing the tube.

drainova® tunneler: Do not use the product if there is visible damage. Sharp edges or splintering could cause irritations or perforation of tissue. The tunneler is only to be used for implanting an indwelling catheter or a comparable implant.

Important: Dispose of used material in accordance with your local, state and national regulations. Where necessary, provide a waste disposal container.

The product may only be implanted by medical professionals! The subsequent drainage can be carried out by persons without medical training.

Please also observe the hygiene regulations of your institution.

Do not use substances that may damage the components.

The catheter must not be cleaned with octenisept® wound disinfectant.

The tearaway introducer must not come into contact with alcohol, acetone or solutions containing these substances.

Ensuring traceability.

For clear and complete traceability, the self-adhesive labels of the catheter set with an unique identification number (device identification) e.g. UDI/LOT must be affixed to the medical prescription, the implant card and the patient file. The self-adhesive labels are located in the lower area of the label on the outer sterile packaging (inside cardboard box).

Possible undesirable side effects can include:

Reexpansion pulmonary oedema*, pneumothorax*, hemothorax, tear or injury of the lung, liver or other organs, hypotension*, circulatory collapse*, wound infection, empyema, seroma, electrolyte imbalance, loss of protein, leakage of fluid from incision site, peritonitis*, infection*, development of metastases, bleeding, septa* in the affected body cavities, chest or abdominal pain, catheter dislocation, hepatorenal syndrome*, catheter blockage, defect catheter, defect valve, complications during explantation (including residual fragments of the catheter) or failed implantation, Skin burns during defibrillation.

The alcohol swabs are highly flammable.

Keep the swabs away from naked flames.

Please note: Effusion drainage can cause some discomfort or be painful for the patient. Should the patient experience severe pain, it is possible to stop or slow down the drainage flow by using the clamp on the connection tube.

Please note: All serious incidents relating to the product must be reported to the manufacturer and relevant authorities.

WARNINGS

regarding implantation

Particular care must be taken to ensure that the catheter remains sterile during implantation.

The implantation of the catheter should be carried out under sterile conditions in an aseptic environment. As soon as the catheter is removed from its sterile packaging, care must be taken to ensure that it does not come into contact with unsterile objects.

The catheter may become electrostatically charged attracting particles from the surrounding environment, which can cling to the silicone. It is therefore important to carry out the implantation in a speedy manner. Correct use minimises risk of infection.

Exercise caution when carrying out the implantation.

When puncturing, take care that no surrounding tissue or organs are injured. Should resistance be experienced at any point during the procedure, the operation should be stopped immediately and the reason for the resistance should be determined using a visual aid.

If after filling the syringe the solution becomes cloudy or discoloration or particles occur, immediately stop use, inform ewimed GmbH and if possible return the syringe to ewimed GmbH. The guidewire (Seldinger wire) must not be pulled back through the cannula as withdrawal can cause damage to the wire. Care must be taken that the guidewire is not inserted deeper than necessary. If resistance is felt during handling with the guide wire, the procedure must be stopped immediately and the reason for the resistance must be determined using various visual examination procedures.

If used improperly, guide wires can lead to perforation of tissue. The consequences of perforation can be severe and can lead to death of the patient.

Exercise caution when handling the scalpel in order to avoid unwanted injuries. When making the incision with the scalpel take care not to damage the guidewire.

Take care not to damage the catheter when placing the sutures. Do not constrict the catheter when placing the suture to hold it in place. Chiraflon Monofilament elicits a minimal initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. It does not interfere with actual wound infection.

It is important to ensure the correct placement of the catheter.

The inserted end of the catheter is perforated. The drainage holes must be completely submerged in the body cavity that requires drainage, otherwise effusion fluid can leak into the subcutaneous tunnel.

The barium sulfate strip may look intermittent in X-rays due to the drainage holes.

EN

WARNINGS

post implantation


EN

Never cut through the catheter or cut the safety valve off.

Ensure that the catheter does not come into contact with any sharp objects.

If the functionality of the safety valve cannot be guaranteed or if the safety valve has been cut off by mistake (i.e. is no longer attached), proceed as follows:

- firmly pinch the catheter closed with your fingers
- using your other hand, hold the emergency clamp open so that the catheter can be fed through the open end of the clamp.
- push the catheter completely through the narrow end of the clamp. Ensure that the emergency clamp is not attached directly on the valve, but at a distance of approx. 2 cm.
- Immediately inform your doctor.

Only use original accessories.

Do not attempt to connect the catheter to non-original drainage products. Only the specially-developed access tip can be inserted into the safety valve. The use of foreign items can damage the safety valve. Improper handling may cause air to enter the body or fluid to leak uncontrollably.

The cone adapter is used solely for connecting compatible drainage systems with open tube, via the drainova® irrigation tube, to the catheter. Ensure correct fit when connecting the luer lock components

Do not pull on the catheter.

Cautionary measures should be taken to ensure that the catheter is not put under any strain, accidentally dislocated or pulled out. Pulling the catheter out can cause bleeding, lead to infection and/or hinder drainage. In these instances, it may be necessary to re-implant the catheter.

Please observe the recommended maximum volume of fluid to be drained.

ewimed does not recommend draining more than 1000 ml pleural effusion or 2000 ml ascites per day. It is in the responsibility of the medical practitioner to recommend larger drainage amounts. The individual condition of the patient should be taken into consideration.

Please note: If the effusion fluid cannot be drained or there is a gradual reduction in drainage output, it is possible that the catheter or the connection tube is blocked. Gently squeeze the catheter and the connection tube manipulating them in a similar way to milking. Should drainage still not be possible, change the drainage reservoir. If the blockage can not be resolved using the methods described above, the catheter can be flushed out with the rinsing tube.

Please note: The drainova® ArgentiC catheter is not transparent. As a result a blockage is not visible.

Always keep the catheter safety valve clean and the access tip on the connected tube sterile.

Take care to ensure that they do not come into contact with any non-sterile objects in order to avoid contamination. The safety valve must be disinfected before each drainage.

While draining, make sure that the safety valve and the access tip are securely connected. Take care not to tug or pull on the connected tube.

Take particular care to ensure that the access tip is inserted in the safety valve in a sterile manner and make sure that the connection is secure. It can lead to contamination if the parts become disconnected by mistake. In case this occurs, disinfect the safety valve with a new alcohol wipe and dispose of the disconnected drainage material. Use a new sterile drainage kit to continue with the drainage.

Please note: If required (for additional protection against accidental disconnection) the connection can be secured with the separately available drainova® clickFix (catalogue number P1100 / P1100S).

The drainova® clickFix can be additionally prescribed. It is not contained in the drainage kit.

Cleaning / disinfection procedure if contaminated.

If effusion fluid is leaking, clean the skin thoroughly with soap and water. Use suitable cleaning or disinfectant products for surfaces. Do not use octenisept® or disinfectants containing iodine for cleaning / disinfecting the catheter or safety valve.

Please note: The catheter is magnetic resonance safe.

EN

GENERAL PREPARATION

These instructions serve as a guideline. The physician carrying out the implantation is responsible for the implementation of the procedures.

The implantation of the catheter should be carried out under local anaesthesia on a sedated patient. Imaging techniques (e.g. ultra sound, colour duplex, X-rays) can be used to aid correct placement of the catheter and to avoid misplaced incisions and injuries to organs and tissue.

The patient's anatomy and general condition play an important role when choosing the catheter implantation site. Ideally, the catheter exit site should be placed so that the patient can easily

access the safety valve, allowing them to carry out drainage and dressing changes independently.

The perforated end of the catheter can be shortened according to the individual anatomy of the patient. It can be cleanly cut off in the space between two holes using for example a scalpel. Care should be taken to ensure the cut is as straight and even as possible.


RECOMMENDED IMPLANTATION PROCEDURE

Pleural Effusion

EN

The most suitable medical and surgical procedure should be determined by a doctor. Which procedure that will be, depends on the individual situation of the patient involved.

The following procedure describes the placement of the catheter in the pleural cavity. The following illustrations show the catheter during implantation in the pleural cavity.

- A** 1. Incision: Puncture site / catheter entrance site
- B** Tunnel
- C** 2. Incision: Catheter exit site
- D** Catheter

1. Check the material is complete prior to starting.
2. Prepare the patient for the procedure. Locate the effusion fluid using ultrasound.

Please note: For a safer implantation with fewer complications it is advantageous to have a large volume of effusion.

3. Mark the incision site (A) and the catheter exit site (C) and determine where the subcutan-

eous tunnel (B) should run. The incision site in the intercostal space is determined by the location of the effusion. The catheter exit site is typically caudal-ventral and marks the end of the approx. 5 cm long tunnel. For better access and patient comfort, the chosen exit site should be positioned ventrally in the same intercostal space.



4. Position the patient for the procedure in a way that the implantation site can be easily accessed.
5. Disinfect a large area for implantation and cover with sterile material.
6. Apply a generous amount of local anesthetic to the puncture incision site (A), the catheter exit site (C) and the tunnel (B).

Please note: The tunnel (B) should be infiltrated with 20 – 30 ml local anaesthetic, 1% fan-shaped. This will lessen the resistance of the subcutaneous tissue* and the tunneler can slide more easily through the tissue. Allow enough time for the anaesthesia to take effect.

7. Carefully push the syringe with attached needle through the chosen intercostal space on the upper edge of the ribs in the pleural gap. Take care not to damage surrounding tissue or organs.
8. If effusion with can be aspirated, leave the needle in position and remove the syringe.
9. Insert the guide wire through the needle as far as necessary. The wire has a j-tip and will not damage surrounding tissue or organs.
10. When the guide wire is in position, retract the needle, whilst holding the guidewire in place.
11. Make a 1,5 cm wide incision on the guidewire parallel to the planned catheter exit site. Take

care not to damage the guidewire with the scalpel.

12. Maintaining a distance of approx 5 cm to the planned catheter exit site, make a second incision (approx 1,5 cm wide) parallel to the first incision.
13. Securely attach the perforated end of the catheter to the intended side of the tunneler.

14. Push the tunneler (E) starting with the tip via the second incision retrogradely, subcutaneously out to the first incision. Pull the catheter (D) through until the polyester cuff is positioned in the middle of the subcutaneous tunnel. Remove the tunneler from the catheter.

Important: Take care that the tunneler passes through the subcutaneous skin tissue and not through the underlying layer of muscular tissue.

Please note: To ensure the correct position of the polyester cuff and to remove any kinks, the catheter will be retracted later on (see point 16). After implantation, the cuff should be positioned approx. 1 cm before the exit from the skin surface in the subcutaneous tunnel. This makes a potential explantation at a later date easier.

15. Insert the 16Fr tearaway introducer over the guidewire in the pleural gap and hold the guidewire in place. Pull out the dialator along with the guidewire whilst holding the 16Fr introducer in position.

Important: In order to prevent effusion leakage make sure to close off the lumen of the introducer with your finger immediately after pulling out the dialator.

Please note: The insertion of the catheter will be impeded if the introducer is bent or kinked.

16. Quickly insert the entire perforated end of the catheter through the introducer in the pleural gap.
17. Break open the introducer using the attached grips and slowly and symmetrically tear it apart. Ensure the catheter is held in place until the introducer has been completely removed.

Please note: During this process, a part of the catheter might be pushed out of the body again. If this happens, simply insert it further in.

EN



38 | Recommended Implantation Procedure

EN

Please note: After removing the introducer there may be a kink in the catheter. This kink can be removed by gently and carefully withdrawing the catheter. Take care that the polyester cuff is not pulled out of the tunnel. The cuff should be positioned approx 1 cm before the catheter exit site. The postoperatively ingrown cuff ensures that the catheter lies securely and helps prevent infection.

Please note: Before applying sutures to the incisions, drain off fluid to check that the implantation was successful. Please refer to chapter 'drainage procedure', page 41.

18. Apply sutures to the first incision at the puncture site (A). Apply sutures to the second incision at the catheter exit site (C) and fix the catheter in place with an additional suture. Take care that the catheter lumen is not constricted.


Please note: The sutures can typically be removed after 7 – 10 days. The suture holding the catheter in place however, should be taken out after approx. 30 days in order to allow the cuff to grow in.

RECOMMENDED IMPLANTATION PROCEDURE

Ascites

Determining the most suitable medical and surgical approach remains the responsibility of the treating physician. Which approach that will be, depends on the individual patient.

The following procedure describes the placement of the catheter in the peritoneal cavity*. The following illustrations show the placing of the catheter in the peritoneal cavity.

- A** 1. Incision: Puncture site / catheter entrance site
- B** Tunnel
- C** 2. Incision: Catheter exit site
- D** Catheter



1. Check the material is complete prior to starting.
2. Prepare the patient for the procedure. Locate the effusion fluid using ultrasound.

Please note: For a safer implantation with fewer complications it is advantageous to have a large volume of effusion.

3. Mark the incision site (A) and the catheter exit site (C) and determine where the subcutaneous tunnel (B) should run. The lateral, abdominal puncture site depends on the location of the ascites. The catheter exit site is typically approx. 8 cm above the puncture site and marks the end of the approx. 5 cm long tunnel. For better access and patient comfort, the chosen exit site should be positioned medially.
4. Position the patient for the procedure in a way that the implantation site can be easily accessed.
5. Disinfect a large area for implantation and cover with sterile material.
6. Apply a generous amount of local anesthetic to the puncture incision site (A), the catheter exit site (C) and the tunnel (B).

Please note: The tunnel (B) should be infiltrated with 20 – 30 ml local anaesthetic, 1% fan-shaped. This will lessen the resistance of the subcutaneous tissue and the tunneler can slide more easily through the tissue. Allow enough time for the anaesthesia to take effect.


7. Carefully push the syringe with attached needle through the chosen peritoneal cavity. Take care not to damage surrounding tissue or organs.
8. If effusion can be aspirated, leave the needle in position and remove the syringe.
9. Insert the guidewire through the needle as far as necessary. The wire has a j-tip and will not damage surrounding tissue or organs.
10. When the guidewire is in position, retract the needle, whilst holding the guidewire in place.
11. Make a 1,5 cm wide incision on the guidewire parallel to the planned catheter exit site. Take care not to damage the guidewire with the scalpel.
12. Maintaining a distance of approx. 8 cm to the planned catheter exit site, make a second incision (approx 1,5 cm wide) parallel to the first incision.
13. Securely attach the perforated end of the catheter to the intended side of the tunneler.

EN

40 | Recommended Implantation Procedure

EN

14. Push the tunneler (E) starting with the tip via the second incision retrogradely, subcutaneously out to the first incision. Pull the catheter through until the polyester cuff is positioned in the middle in the subcutaneous tunnel. Remove the tunneler from the catheter.

Please note: Take care that the tunneler passes through the subcutaneous skin tissue and not through the underlying layer of muscular tissue.

Please note: To ensure the correct position of the polyester cuff and to remove any kinks, the catheter will be retracted later on (see point 16). After implantation, the cuff should be positioned approx. 1 cm before the exit from the skin surface in the subcutaneous tunnel. This makes a potential explantation at a later date easier.

15. Insert the 16Fr tearaway introducer over the guidewire in the peritoneal gap and hold the guidewire in place. Pull out the dialator along with the guidewire, whilst holding the 16Fr introducer in position.

Please note: In order to prevent effusion leakage, make sure to close the lumen of the introducer with your finger immediately after pulling out the dialator.

Please note: The insertion of the catheter will be impeded if the introducer is bent or has kinks in it.


16. Quickly insert the entire perforated end of the catheter through the introducer in the peritoneal cavity.
17. Break open the introducer using the attached grips and slowly and symmetrically tear it apart. Ensure the catheter is held in place until the introducer has been completely removed.

Please note: During this process, a part of the catheter might be pushed out of the body again. If this happens, simply insert it further in.

Please note: Removing the introducer may cause a kink in the catheter. This kink can be removed by gently and carefully pulling the catheter. Take care that the polyester cuff is not pulled out of the subcutaneous tunnel. The cuff should be positioned approx. 1 cm before the catheter exit site. The postoperatively ingrown cuff ensures that the catheter lies securely and helps prevent infection.

Important: Before applying sutures to the incisions, drain off fluid to check that the implantation was successful. Please refer to chapter 'drainage procedure', page 41.

18. Apply sutures to the first incision at the puncture site (A). Apply sutures to the second incision at the catheter exit site (C) and secure the catheter in place with an additional suture. Take care that whilst so doing the catheter lumen is not constricted.

Please note: The sutures can typically be removed after 7 – 10 days. The suture securing the catheter in position should however, be left in place for approx. 30 days in order to allow the cuff to grow in.

EN

DRAINAGE PROCEDURE

The following ewimed clinic drainage products are available for use with the catheter:

- 2000 drainova® reservoir, soft vacuum
- 50-7220 ewimed gravity reservoir - 2000ml
- 50-7210 PleurX™ cacuum reservoir
- P8531 ewimed Pneu-Pack I

Drainage Consumables

- 6060 flushing line
- P1100 drainova® clickFix
- 9060 drainova® silicone cap, sterile

The following ewimed drainage products are available for home care:

- 2010 drainova® drainova® reservoir, drainage set
- 50-7505 ewimed drainage kit, 2000 ml
- 50-7510 PleurX™ drainage kit
- 50-7500B PleurX™ drainage kit


42 | Completion of Drainage Procedure

EN

Important: When using a suction device, set a vacuum capacity of max. -60 cm H₂O or a max. flow velocity of 400 ml per minute. Under certain circumstances the catheter can become stuck. If this is painful for the patient, close the clamp on the connected tube. In so doing the flow can be slowed or completely stopped.

Please observe the corresponding instruction manuals when draining with the above-named products.


All products named above feature an access tip on the connected tube to allow easy and secure connection with the safety valve. A click can be heard and felt when connecting the access tip.



Using the drainova® clickFix, it is possible to secure the connection between the access tip and safety valve. The connection is inserted so that it lies in the specially designed inner ridges and secured in place. Accidental disconnection or contamination can be avoided when the drainova® clickFix is closed.

COMPLETION OF DRAINAGE PROCEDURE

1. Once drainage has been completed, place the drainova® silicone , sterile cap on the safety valve.
2. Place the foam pad with the slit facing up around the catheter, loop the catheter and position it on the foam pad.
3. Cover the catheter with the gauze pad and apply the self-adhesive dressing.

2. Place the foam pad with the slit facing up around the catheter, loop the catheter and position it on the foam pad.

CONNECTING THE RINSING TUBE AND FLUSHING THE CATHETER (DISLODGING OCCLUSIONS)

The following material is required to flush the catheter:


- physiological saline solution (0,9% NaCl)
- Three-way valve
- Syringe
- Cannula 1G
- Alcohol wipe

1. Open a new package of rinsing tubes. The rinsing tube in the package is sterile. It is important that the access tip on the drainova® Rinsing Tube remains sterile.
2. Connect a three-way valve onto the Luer-lock connection of the rinsing tube.
3. Fill a 10ml or 20ml syringe (cannula 1G) with 10 – 20ml physiological saline solution.

4. Connect the syringe onto the three-way valve.
5. Rotate the three-way valve so that the connection between the syringe and the rinsing tube is open.
6. Remove any air bubbles from the rinsing tube. To do this, hold up the end of the rinsing tube and carefully inject the saline solution into the tube until solution runs out.
7. Firmly hold the access tip on the grip. By gently rotating, carefully remove the protective cover from the access tip and discard.

8. Insert the access tip in the safety valve. Take care that the safety valve and the access tip are securely connected when draining. A click can be heard and felt when connected.

9. Flush the catheter by applying pressure to the plunger of the syringe.
10. If necessary, draw up more saline solution with the syringe.
 - 10.1 Squeeze the clamp on the rinsing tube closed.

- 10.2 Close the three-way valve
- 10.3 Remove the syringe from the three-way valve, fill with saline solution and connect it to the three-way valve.
- 10.4 Rotate the three-way valve to the open position and carefully inject saline solution until fluid runs out of the three-way valve.
- 10.5 Adjust the three-way valve as described in step 5 and continue from point 9

11. When the occlusion has successfully been dislodged, complete the flushing procedure by squeezing the clamp closed on the rinsing tube.
12. Afterwards the physiological saline solution should be aspirated and the effusion drained with a product listed on page 41.

Please note: The flushing procedure was successful where both a smooth application of saline solution as well as aspiration of effusion was possible.

EN

CATHETER EXPLANTATION

Possible indications for the removal of the catheter:

- Pleurodesis/no ascites
- Pleural empyema / Peritonitis with strong viscous/ putrid effusion (if treatment with antibiotic therapy was not effective)
- Infection in the subcutaneous tunnel, persisting despite systematic antibiotic therapy
- The catheter is blocked and the occlusion can not be flushed
- The catheter is in a cavity that no longer fills with effusion

- No re-accumulation of effusion fluid determined during preceding sonographic monitoring and drainage-free interval (no specification, typical no-drainage period of 1–4 weeks where patient shows no symptoms).
- The catheter has become displaced

Important: The indication should be checked again at the end of the drainage-free period via ultrasound before removing the catheter.

44 | Catheter Explantation

Example Procedure for Catheter removal:

1. Check the coagulation parameters
2. Prepare a sterile workspace
3. Position the patient to allow easy access to the catheter
4. Infiltrate the skin around the catheter exit site with a generous amount of local anesthetic (observe application time)
5. Externally feel for the polyester cuff (bulge-like thickening)
6. Incision in the skin directly at the catheter exit site
7. blunt, circular freeing of the polyester cuff using the irrigation cannula, dissecting scissors or similar instrument, (the catheter will adhere to the subcutaneous tissue only at the cuff, difficult conditions can arise when removing the catheter either due to the position itself or the walls of the tumor)
8. Use one hand to carefully apply pressure to the subcutaneous tunnel and use the other hand to pull out the catheter. If the catheter as a whole cannot be mobilised it might be necessary to make an additional incision (corresponding to the length, at the end of the tunnel/ using the incision at the puncture site as a guide) in order to recover the part of the catheter that is buried.
9. Subsequent closure of the incision(s) with skin sutures
10. Apply bandage.

EN

GLOSSARY

Ascites	Effusion or accumulation of fluid in the abdomen (stomach)
Circulatory collapse	Temporary loss of consciousness
DEHP	Plasticizers in plastics
Electrolyte imbalance	Disturbance of the balance of ions (e.g. sodium, potassium, calcium) in the blood or cells in the body
Empyema	Infection of the pleura
Hepatorenal syndrome	Decrease in renal function
Hypotension	Low blood pressure
Infection	Reaction of the body's own defence system to pathogens such as bacteria, viruses etc.
Peritoneal cavity	Abdominal cavity
Protein loss	Reduction of the amount of protein in the blood (protein depletion)
Peritonitis	Inflammation of the peritoneum (abdominal cavity)
Phthalate	Plasticizer (e.g. DEHP)
Pleural effusion	Effusion or fluid accumulation in the pleural cavity
Pleurodesis	Adhesion of the pleural space
Pneumothorax	Air in the space between the lungs and the chest wall
Reexpansion pulmonary edema	Accumulation of fluid in functional lung tissue due to the sudden expansion of the previously compressed lung; acute shortness of breath
Septa	Tissue bridges / chambering between separated spaces in a body cavity (e.g. abdomen)
Subcutaneous tissue	Tissue under the skin

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PleurX™/PeritX™ is a trademark and/or registered trademark of Becton Dickinson and Company, or one of its subsidiaries.
drainova® is a trademark and/or registered trademark of ewimed GmbH.
octenisept® is a trademark of Schülke & Mayr GmbH.
Germany | ewimed GmbH | Dorfstrasse 24 | 72379 Hechingen | Tel.: +49 7471 73972-0 | info@ewimed.com