DE drainova® Katheter-Set / drainova® ArgentiC (S) Katheter-Set / **EN** drainova® catheter set / drainova® ArgentiC (S) catheter set / **BG** Комплект катетър drainova® / комплект катетър drainova® ArgentiC (S) / **DA** drainova® katetersæt / drainova® ArgentiC (S) katetersæt / **EL** σετ καθετήρων drainova® / σετ καθετήρων drainova® ArgentiC (S) / **ES** Juego de catéter drainova® / juego de catéter drainova® ArgentiC (S) / **ET** drainova® kateetri komplekt / drainova® ArgentiC (S) kateetri komplekt / **FI** drainova® katetrisarja / drainova® ArgentiC (S) katetrisarja (S) / **FR** Kit de cathéter drainova® / Kit de cathéter drainova® ArgentiC (S) / **HR** Komplet katetera drainova® / komplet katetera drainova® ArgentiC (S) / **HU** drainova® katéterkészlet/drainova® ArgentiC (S) katéterkészlet / **IT** drainova® Set di cateteri / drainova® ArgentiC Set di cateteri (S) / **LT** drainova® kateterių rinkinys / drainova® ArgentiC (S) kateterių rinkinys / **LV** drainova® katetru komplekts / drainova® ArgentiC (S) katetru komplekts / **MT** Sett tal-kateter drainova® / Sett tal-kateter drainova® ArgentiC (S) / **NL** drainova® katheterset/drainova® ArgentiC (S) katheterset / **NO** drainova® katetersett / drainova® ArgentiC (S)-katetersett / **PL** Zestaw cewników drainova® / zestaw cewników drainova® ArgentiC (S) / **PT** Conjunto de cateteres drainova® / Conjunto de cateteres drainova® ArgentiC (S) / **RO** Set pentru cateter drainova®/set pentru cateter drainova® ArgentiC (S) / **SK** Súprava katétra drainova®/súprava katétra drainova® ArgentiC (S) / **SL** Komplet katetra drainova®/komplet katetra drainova® ArgentiC (S) / **SV** drainova® kateterset / drainova® ArgentiC (S) kateterset / **TR** drainova® Kateter Seti / drainova® ArgentiC (S) Kateter Seti

DE Gebrauchsanweisung / **EN** Instructions for use / **BG** Инструкции за употреба / **DA** Brugsanvisning / **EL** Οδηγίες χρήσης / **ES** Instrucciones de uso / **ET** Kasutusjuhend / **FI** Käyttöohjeet / **FR** Mode d'emploi / **HR** Upute za uporabu / **HU** Használati utasítás / **IT** Istruzioni per l'uso / **LT** Naudojimo instrukcijos / **LV** Lietošanas instrukcija / **MT** Struzzjoni għall-użu / **NL** Gebruiksaanwijzing / **NO** Bruksanvisning / **PL** Instrukcja obsługi / **PT** Instruções de utilização / **RO** Instrucţiuni de utilizare / **SK** Návod na použitie / **SL** Navodila za uporabo / **SV** Instruktioner för användning / **TR** Kullanım Kılavuzu

2 DE KONTAKT / EN CONTACT US / BG КОНТАКТ / DA KONTAKT / EL ΣΤΟΙΧΕΙΑ ΕΠΙΚΟΙΝΩΝΙΑΣ / ES CONTACTO / ET KONTAKTANDMED / FI YHTEYS / FR CONTACT / HR KONTAKT / HU KAPCSOLAT / IT CONTATTI / LT KONTAKTINĖ INFORMACIJA / LV KONTAKTINFORMĀCIJA / MT KUNTATT / NL CONTACT / NO KONTAKT / PL KONTAKT / PT CONTACTO / RO CONTACT / SK KONTAKT / SL KONTAKTNI PODATKI / SV KONTAKTA OSS / TR İLETİŞİM

AUSTRIA

ewimed austria GmbH
Tullnerbachstrasse 92a
3011 Neu-Purkersdorf
Tel.: +43 2231 2250-0
info@ewimed.com

BELGIUM

ewimed Belgium BV
Kortrijksesteenweg 1126a
9051 Sint-Denijs-Westrem
Tel.: +32 7848 0510
Info.BE@ewimed.com

DENMARK

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ENGLISH

Please read the instructions for use carefully before using the drainova® catheter set / drainova® ArgentiC (S) catheter set.

All explanations of terms (*) can be found at "Glossary" on page 24.

THE DRAINOVA® CATHETER SET / DRAINOVA® ARGENTIC (S) CATHETER SET

- Used for the implantation of a drainova® catheter / drainova® ArgentiC (S) catheter for recurrent, refractory, malignant, and non-malignant effusion collections in serous body cavities.
- Relief of symptoms through drainage of accumulated fluids.

Indication

The use of the catheter is indicated in the following situations:

- Malignant pleural effusion*
- Non-malignant pleural effusion*
- Malignant ascites*
- Non-malignant ascites*

20 Contraindication

The use of the catheter is contraindicated in the following situations:

- The body cavity has multiple chambers (septa*), so that no relief of dyspnea or other symptoms can be expected after drainage.
- Coagulopathy
- Infection* of the body cavity
- Chylous effusion
- Displacement of the mediastinum by more than 2 cm to the ipsilateral side of the pleural effusion*.
- Allergy to product materials.

The use of catheter set components is contraindicated in the following situations:

- Syringe: use only for the intended purpose, not suitable for highly viscous fluids.
- Guide wire contains cobalt: Since the guide wires are used as a medical device for the applications mentioned above, the contraindications for these applications also essentially apply:
 - Central nervous system applications.
 - Applications on the central circulatory system.
- Introducer: Increased risk of pneumothorax in patients with severe chronic lung disease. Deteriorated wound healing possible with previously irradiated anterior thoracic wall.

INTENDED USE

The drainova® catheter set / drainova® ArgentiC catheter set serves as a comprehensive set of components required for implantation of a drainova® catheter / drainova® ArgentiC catheter in recurrent, refractory, malignant, and non-malignant effusion accumulations in serous body cavities to relieve the symptoms by draining the fluids.

This includes:

- Malignant pleural effusion*
- Non-malignant pleural effusion*
- Malignant ascites*
- Non-malignant ascites*

All set components are CE-marked and sterilized with ethylene oxide gas.

THE DRAINOVA® CATHETER

- Material: The catheter is made of biocompatible silicone. **1**
- Length and diameter: 66 cm long, diameter 15.5 Fr.
- Perforated area: 26 cm with opposite holes for fluid inflow (implanted part).
- Safety valve at the non-implanted end of the catheter.
 - Closes automatically when disconnected. Enables connection with drainage accessories via access tip.
 - Prevents air ingress and fluid leakage, as well as the unwanted escape of effusion.
 - May only be opened with suitable materials that have the specially developed access tip.
- Polyester cuff between the perforated side and the safety valve.
 - Fixes the catheter in the subcutaneous tunnel with the skin tissue.
- Radiopaque strip
 - Barium sulphate strip over the entire length of the catheter for X-ray visibility.

THE DRAINOVA® ARGENTIC (S) CATHETER

- Geometry and handling: Identical to the drainova® catheter. **1**
- Material: Biocompatible silicone with embedded microsilver.
- Function:
 - Silver protects against bacterial colonization on surfaces.
 - Reduces the risk of infection*.

The drainova® ArgentiC S catheter is a variant of the proven drainova® ArgentiC catheter and is characterized by a shortened tube length. **1**

WARNINGS

General

Do not use the product if:

- the sterile packaging is damaged
- the product is damaged
- the expiry date has passed

The product is to be used exclusively for its intended purpose.

Patient target group: The medical device may only be used in adults.

Intended use environment: The catheter set is intended for use in clinical facilities (such as hospitals and medical practices).

The consumables are intended for single use only!

- Reuse may contribute to cross-contamination.
- The catheter is intended for one-time implantation.
- The catheter set is delivered in a sterile condition, which is ensured by suitable quality controls.
- The catheter sets must not be resterilized.
- ewimed GmbH accepts no liability for products that have been reused, resterilized, or removed from damaged packaging.

drainova® clickFix **2**: Insert the tube connection correctly into the drainova® clickFix to prevent the tube from being crushed.

drainova® tunneler: The tunneler contains 8 - 10 % nickel.

- Do not use if there is visible damage.
- Sharp edges / chipping can cause irritation or tissue perforation.
- The tunneler is to be used exclusively for tunneling a long-term catheter to be implanted or a comparable implant.
- The tunneler is supplied sterile.
- Do not use it non-sterile, as this can cause infections.

Note: When disposing of used material, please observe the local, state, and national regulations. If necessary, provide a residual waste container.

User group: The product may only be implanted by qualified medical professionals! The subsequent drainage may be carried out by persons without a medical background.

Also observe the hygiene regulations of your facility.

Do not use any substances that could damage the components.

- Do not clean the catheter / safety valve with octenisept® or a disinfectant containing iodine.
- The introducer must not come into contact with alcohol, acetone, or solutions containing these substances.
- Step adapter:
 - Only use alcohol-based disinfectants without chemical additives.
 - Only use disinfectant when necessary.
 - After cleaning, ensure that the connected products are correctly and firmly seated.

Ensure traceability:

- For clear and complete traceability, the self-adhesive catheter set labels with unique identification number (product identification), e.g. UDI / LOT, must be affixed to the doctor's prescription, the implant card, and the patient file.
- The self-adhesive labels are located at the bottom of the label on the outer sterile packaging (inside the cardboard box).

Possible undesirable side effects may include:

Re-expansion pulmonary edema*, pneumothorax*, hemothorax*, laceration/injury to lung, liver, or other organs, hypotension*, circulatory collapse*, wound infection, wound healing disorders, traumatized tissue, empyema*, seroma, electrolyte imbalance*, protein loss*, leakage of body fluid at the incision, peritonitis*, infection*, spread of metastases, bleeding, septa* in the affected body cavities, pain in the chest/abdominal area, catheter dislocation, hepatorenal syndrome*, catheter blockage, catheter malfunction, valve malfunction, complication during explantation (esp. catheter fragments left behind), or failed implantation, skin burns during defibrillation.

The alcohol swabs are highly flammable.

The swabs must be kept away from naked flames.

Note: Draining the effusion can be unpleasant or even painful for the patient. If the pain is too severe, the drainage flow can be stopped or slowed down for a few minutes by pressing the clamp on the connecting tube.

Note: All serious incidents in connection with the product must be reported to the manufacturer and the responsible authority. Observe the local disposal regulations.

WARNINGS**For implantation**

Care must be taken to ensure that the catheter remains sterile during implantation.

- Implantation: Insert the catheter in an aseptic room in compliance with sterile precautions.
- Catheter: Do not allow contact with non-sterile objects
 - This can attract airborne particles that adhere to the silicone (electrostatic charge).
 - Rapid placement of the catheter is required.
- Proper use minimizes the risk of infection.

Take particular care during implantation.

- Surrounding tissue and organs must not be damaged.
- If the cannula shows resistance, stop the procedure and visually check the cause.
- After filling the syringe, check for changes in the solution:
 - If the solution becomes cloudy, discolored, or contains particles, discontinue use.
 - Inform ewimed GmbH and return the syringe if necessary.
- Do not pull back the guide wire (Seldinger wire), as this can damage the wire.
 - Only insert the wire as deep as necessary.
 - If resistance is felt, stop the procedure and clarify the cause with imaging procedures.
 - If used incorrectly, the guide wire can lead to tissue perforation. The consequences of perforation can be serious and lead to death.
- Handle the scalpel with care to avoid cuts.
- Do not damage the guide wire during the incision.
- Do not damage or constrict the catheter during skin suture and anchoring suture.
- Chirafilon monofilament causes a minimal inflammatory reaction with subsequent fibrotic encapsulation. No impairment of wound healing.

Care must be taken to ensure correct placement of the catheter.

- The inserted end of the catheter is perforated.
- The drainage holes must be completely in the body cavity to be drained. Otherwise, effusion may enter the subcutaneous* tunnel.
- The barium sulphate strip may appear interrupted by the drainage holes in the X-ray image.

WARNINGS**After implantation**

Never cut through the catheter or cut off the safety valve.

Avoid contact with sharp objects.

If the function of the safety valve is no longer guaranteed or if the safety valve has been accidentally cut off or is no longer present, proceed as follows: 3

- Press the catheter firmly closed with your fingers.
- Open the clamp with one hand and insert the catheter into the open end.
- Use a clamp available at the clinic. Position the clamp at a distance of approx. 2 cm from the valve.
- Contact the attending physician immediately.

Only use compatible products.

- Only compatible products may be connected to the catheter.
- Nothing other than the access tip specially developed for this purpose may be inserted into the safety valve; otherwise the safety valve may be damaged.
- Improper handling can lead to air entering the body or uncontrolled discharge of effusion.
- The step adapter is used exclusively to connect suitable drainage systems with an open tube to the catheter via the drainova® flushing line.
- When connecting the Luer-Lock components, ensure that they are correctly seated.

Do not pull on the catheter.

- Precautions must be taken to ensure that the catheter is not pulled, accidentally dislocated, or pulled out.
- Pulling out can lead to bleeding and infection* and/or prevent drainage.
- It may be necessary to re-implant the catheter.

Observe the maximum recommended drainage volume.

- ewimed GmbH recommends draining no more than 1000 ml pleural effusion* or 2000 ml ascites* per day.
- Draining larger amounts of effusion is the responsibility of the attending physician.
- Consider the individual condition of the patient.

Note: If no more effusion is drained or the amount gradually decreases, the catheter or connection tube may be blocked. Gently squeeze the catheter and the connecting tube, as if milking. If the drainage does not start, change the inserted reservoir. If a catheter blockage cannot be cleared with the above measures, it can be flushed with the flushing line.

Note: The drainova® ArgentiC (S) catheter is not optically transparent. This means that a blockage in the catheter cannot be seen.

Keep the safety valve on the catheter clean and the access tip on the connection tube sterile.

- Do not allow contact with non-sterile objects, in order to avoid soiling or contamination.
- The safety valve must be disinfected before each drainage.

Ensure that the safety valve and the access tip are fully connected during the drainage process. Ensure that the connection tube is not tugged or pulled.

- Ensure a sterile procedure when connecting the access tip to the safety valve.
- The plug connection must be firmly connected.
- Risk of contamination during disconnection. In such a case, disinfect the safety valve with a fresh alcohol swab.
- Dispose of disconnected drainage material.
- A new, sterile drainage set must be used for continuation.
- Observe the local disposal regulations.

Note: If required, the plug connection can be secured with a separately available drainova® clickFix (catalog no. P1100 / P1100S) for additional protection against accidental disconnection.

The drainova® clickFix can be additionally prescribed and is not part of the catheter set. **2**

Cleaning/disinfection in case of contamination.

- In case of effusion leakage: Clean skin with soap and water.
- For surfaces: Use a suitable cleaning agent or disinfectant.
- Catheter / safety valve: Do not use octenisept® or disinfectant containing iodine.

Note: The catheter is magnetic resonance safe.

GENERAL PREPARATORY MEASURES

These instructions are intended as a guide; the implementation of the procedures is the responsibility of the implanting physician.

- Implantation:
 - Perform under sedation and local anesthesia.
 - Use imaging techniques (e.g. ultrasound, color duplex, fluoroscopy) to determine the location of the catheter and avoid misplaced punctures or organ/vascular injuries.
- Selection of the implantation site:
 - Consider patient anatomy and patient condition.
 - Select the exit site so that the patient has access to the catheter safety valve and can independently perform intermittent drainage and dressing changes.
- Adjustment of the catheter:
 - The perforated end of the catheter can be shortened individually. For example, use a sterile scalpel to cut it straight and as flat as possible between two openings.

RECOMMENDED IMPLANTATION PROCEDURE

Pleura*

The appropriate medical and surgical procedure is the responsibility of the doctor. The appropriate procedure depends on the individual patient's situation.

The following procedure describes the placement of the catheter in the pleural cavity. **4**

1. Before starting, check the materials available for completeness
2. Prepare the patient for the procedure. Use ultrasound to localize the position of the effusion.

Note: A large effusion is advantageous for safe implantation with few complications.

3. Select the puncture site (A) and catheter exit site (C) and determine the tunnel route (B). The puncture site in the intercostal space depends on the localization of the effusion. The exit site of the catheter is usually located caudal-ventral and marks the end of the approx. 5 cm tunnel section. For better accessibility and patient comfort, the exit site should be chosen in the ventral direction in the same intercostal space.
4. Position the patient for the procedure so that the implantation site is freely accessible.
5. Disinfect the implantation area extensively and cover it with a sterile dressing.
6. Generously anesthetize the puncture site (A), catheter exit site (C), and tunnel route (B) locally.

Note: The tunnel route (B) should be infiltrated with 20 - 30 ml of local anesthetic, 1 % in a fan-shaped infiltration. This reduces the resistance of the

subcutis, and the tunneler can glide through the tissue more easily. Observe the application time.

7. Advance the cannula with the syringe attached through the selected intercostal space at the upper edge of the rib into the pleural cavity while aspirating. Surrounding tissue and organs must not be damaged.
8. If effusion can be aspirated, leave the cannula in position 15 and withdraw the syringe.
9. Insert the guide wire through the cannula as far as necessary. The wire has a J-tip and does not damage surrounding tissue and organs.
10. Once the guide wire is in position, the cannula is withdrawn. Hold the guide wire in the desired position.
11. Make an approx. 1.5 cm wide incision on the guide wire parallel to the planned exit site of the catheter. Do not damage the guide wire with the scalpel.
12. Make a second incision approx. 1.5 cm wide parallel to the first incision at a distance of approx. 5 cm from the planned exit site of the catheter.
13. Attach the perforated end of the catheter firmly to the side of the tunneler intended for this purpose. **5**
14. Push the tunneler (E) with the tip retrograde via the second incision, subcutaneously out towards the first incision. Retract the catheter (D) until the polyester cuff is positioned in the center of the subcutaneous tunnel. Remove the tunneler from the catheter. **6**

Note: Ensure that the tunneler is passed through the subcutaneous tissue* and not through the underlying muscular layer.

Note: The catheter is later retracted to ensure the correct position of the polyester cuff and to compensate for any kinking of the catheter (see point 16). After implantation, the cuff should be located approx. **1 cm** before it emerges from the skin level in the subcutaneous tunnel. This also facilitates any subsequent explantation.

15. Insert the 16 Fr introducer over the guide wire into the pleural cavity. Pull out the guide wire together with the stylet, leaving the 16 Fr introducer sheath in position.

Note: Immediately after pulling out the stylet, close the lumen of the introducer sheath with your finger to prevent effusion from escaping.

Note: The introducer sheath must not be bent or kinked, as this makes it more difficult to insert the catheter.

16. Push the perforated end of the catheter through the introducer sheath quickly and completely into the pleural cavity.
17. Break open the introducer sheath using the handles provided and pull it apart symmetrically and slowly. Hold the catheter in position until the introducer sheath is completely pulled out.

Note: The catheter can be pushed out of the body again a little; in this case, push the catheter further forward.

Note: Once the introducer sheath has been removed, the catheter may kink in this area. This kink can be removed by carefully withdrawing the catheter. Make sure that the polyester cuff is not pulled out of the tunnel. The polyester cuff should be placed at the end approx. **1 cm** before the exit point of the catheter (C). The cuff, which grows in during the postoperative period, guarantees a secure catheter placement and acts as an infection barrier.

Note: Before suturing the incisions, use drainage to check the success of the implantation. See the chapter "Drainage procedure" on page 23

18. Suture the first incision at the insertion site (A) with a skin suture. Suture the second incision at the exit site (C) and fix the catheter (D) with an anchoring suture. Ensure that the lumen of the catheter is not constricted. **7**

Note: The skin suture can normally be removed after 7 - 10 days. The retaining thread (anchoring suture) should normally be left in place for 30 days to ensure that the cuff grows in.

RECOMMENDED IMPLANTATION PROCEDURE

Ascites*

The appropriate medical and surgical procedure is the responsibility of the doctor. The appropriate procedure depends on the individual patient's situation.

The following procedure describes the placement of the catheter in the pleural cavity. **8**

1. Before starting, check the materials available for completeness
2. Prepare the patient for the procedure. Use ultrasound to localize the position of the effusion.

Note: A large effusion is advantageous for safe implantation with few complications.

3. Select the puncture site (A) and catheter exit site (C) and determine the tunnel route (B). The lateral, abdominal puncture site depends on the localization of the effusion. The exit site of the catheter is usually approx. 8 cm above the puncture site and marks the end of the tunnel section. For better accessibility and patient comfort, the exit site should be chosen in a medial direction. **8**
4. Position the patient for the procedure so that the implantation site is freely accessible.
5. Disinfect the implantation area extensively and cover it with a sterile dressing.
6. Generously anesthetize the puncture site (A), catheter exit site (C), and tunnel route (B) locally.

Note: The tunnel route (B) should be infiltrated with 20 - 30 ml of local anesthetic, 1% in a fan-shaped infiltration. This reduces the resistance of the subcutis, and the tunneler can glide through the tissue more easily. Observe the application time.

7. Carefully advance the cannula with the syringe attached into the peritoneal cavity* while aspirating. Surrounding tissue and organs must not be damaged.
8. If effusion can be aspirated, leave the cannula in position and remove the syringe.
9. Insert the guide wire through the cannula as far as necessary. The wire has a J-tip and does not damage surrounding tissue and organs.
10. Once the guide wire is in position, withdraw the cannula. Hold the guide wire in the desired position.
11. Make an approx. 1.5 cm wide incision on the guide wire parallel to the planned exit site of the catheter (C). Do not damage the guide wire with the scalpel.
12. Make a second, approx. 1.5 cm wide incision parallel to the first incision at a distance of approx. 8 cm from the planned exit point of the catheter (C).
13. Attach the perforated end of the catheter firmly to the side of the tunneler intended for this purpose. **5 9**
14. Push the tunneler (E) with the tip via the second incision retrograde, subcutaneously out towards the first incision. Retract the catheter (D) until the polyester cuff is positioned in the center of the subcutaneous tunnel. Remove the tunneler from the catheter **10**

Note: Ensure that the tunneler is passed through the subcutaneous tissue* and not through the underlying muscular layer.

Note: The catheter is retracted later to ensure correct positioning of the polyester cuff and to compensate for any kinking of the catheter. (See point 16). After implantation, the cuff should be located approx. **1 cm** before it emerges from the skin level in the subcutaneous tunnel. This also facilitates any subsequent explanation.

15. Insert the 16 Fr introducer through the guide wire into the peritoneal cavity* and hold the guide wire in position. Pull out the stylet together with the guide wire, leaving the 16 Fr introducer sheath in position.

Note: Immediately after pulling out the stylet, close the lumen of the introducer sheath with your finger to prevent effusion from escaping.

Note: The introducer sheath must not be bent or kinked, as this makes it more difficult to insert the catheter.

16. Quickly push the catheter with the perforated end through the introducer sheath completely into the peritoneal cavity*.
17. Break open the introducer sheath using the handles provided and pull it apart symmetrically and slowly. Hold the catheter in position until the introducer sheath is completely pulled out.

Note: The catheter can be pushed out of the body again a little; in this case, push the catheter further forward.

Note: Once the introducer sheath has been removed, the catheter may kink in this area. This kink can be removed by carefully withdrawing the catheter. Make sure that the polyester cuff is not pulled out of the tunnel. The polyester cuff should be placed at the end approx. **1 cm** before the exit point of the catheter (C). The cuff, which grows in during the postoperative period, guarantees a secure catheter placement and acts as an infection barrier.

Note: Before suturing the incisions, use drainage to check the success of the implantation. See the chapter on drainage procedures.

18. Suture the first incision at the insertion site (A) with a skin suture. Suture the second incision at the exit point of the catheter (C) and fix the catheter in place using an anchoring suture. Ensure that the lumen of the catheter is not constricted. **10**

Note: The skin suture can normally be removed after 7 - 10 days. The retaining thread (anchoring suture) should normally be left in place for 30 days to ensure that the cuff grows in.

DRAINAGE PROCEDURE

The following ewimed products are available to the clinic for drainage with the catheter:

- 2000 drainova® reservoir, soft vacuum
- 50-7220 ewimed gravity reservoir - filling capacity 2000 ml
- 50-7210 PleurX™ vacuum-tight reservoir
- P8531 ewimed Pneu-Pack I
- 6060 drainova® flushing line
- P1100 drainova® clickFix
- 9060 drainova® silicone cap, sterile

The following ewimed products are available for drainage in home care:

- 2010 drainova® reservoir, drainage set
- 50-7505 ewimed drainage set, 2000ml
- 50-7510 PleurX™ drainage set
- 50-7500B PleurX™ drainage set

Note: When a suction device is used, a maximum suction capacity of -60 cm H₂O or a maximum flow rate of 400 ml per minute must be set. The catheter may get stuck under certain circumstances. If this is painful for the patient, close the clamp on the connection tube. This can slow down or stop the flow of fluid.

When draining with the above-mentioned products, please observe the respective instructions for use.

All the products listed above have an access tip on the connection tube, which can be used to connect the safety valve easily and securely. When the access tip engages, a characteristic click can be heard and felt. **11**

- The connection between the access tip and the safety valve can be fixed with the drainova® clickFix.
- Insert the connection into the inner ribs of the drainova® clickFix. **2**
- The connection is fixed by closing the drainova® clickFix.
 - Prevents disconnection of the connection.
 - Protects the plug connection from soiling.

ENDING THE DRAINAGE PROCESS

1. After successful drainage, place the drainova® silicone cap, sterile, on the safety valve. **12**
2. Place the foam slit compress around the catheter with the slit facing upwards, roll up the catheter, and place it in the foam slit compress. **13, 14**
3. Cover the catheter with the gauze compresses and apply the self-adhesive film dressing. **15, 16**

CONNECTING THE FLUSHING LINE AND FLUSHING THE CATHETER (FLUSHING FREE OF OCCLUSIONS)

The following materials are also required for an irrigation procedure:

- Physiological saline solution (0.9 % NaCl)
 - Three-way stopcock
 - Syringe
 - Cannula 1G
 - Alcohol swab
1. Open the packaging of the flushing line. The drainova® rinsing tubing is sterile, and care must be taken to ensure that it remains sterile.
 2. Connect the three-way stopcock to the Luer-Lock connection of the flushing line.
 3. Fill a 10 ml or 20 ml syringe (1 G cannula) with 10 - 20 ml of a physiological saline solution.
 4. Connect the syringe to the three-way stopcock.
 5. Adjust the three-way stopcock so that the connection between the syringe and the irrigation tube is open.
 6. Bleed the flushing line until there are no bubbles. Lift the end of the tube and fill in saline solution until liquid emerges.
 7. Grasp the access tip at the holding area, carefully pull off the protective cover, and dispose of it. **17**
 8. Push the access tip into the safety valve and ensure that both are firmly connected. You will hear and feel it click into place. **11**
 9. Flush the catheter by pressing on the plunger of the syringe.
 10. If necessary, draw up saline solution again with the syringe.
 - 10.1 Close the squeeze clamp of the flushing line by pressing it together.
 - 10.2 Close the three-way stopcock.
 - 10.3 Disconnect the syringe from the three-way stopcock, fill it with saline solution, and reconnect it.
 - 10.4 Open the three-way stopcock and carefully inject saline solution until liquid emerges.
 - 10.5 Re-adjust the three-way stopcock as described in step 5 and continue with step 9.
 11. To complete the rinsing process, after successfully removing the occlusion, close the pinch clamp on the flushing line by squeezing it.
 12. The saline solution should then be aspirated again and the effusion drained with a product listed under "Drainage procedure" on page 23.

Note: If the saline solution can be applied easily and the effusion can be aspirated easily, the irrigation procedure was successful.

EXPLANATION OF THE CATHETER

Possible indications for removal of the catheter:

- Pleurodesis* / no ascites*
- Pleural empyema / peritonitis* with highly viscous / putrid effusion (if treatment with antibiotics has failed)
- Infection* of the tunnel route that persists despite systemic antibiotics
- Catheter is blocked and the occlusion cannot be cleared
- Catheter is in a chamber that no longer fills with effusion
- No effusion formation with previous sonographic control and drainage-free interval (no specifications, usually drainage break of 1 - 4 weeks if the patient is asymptomatic)
- Catheter is dislocated

Note: Before the catheter is removed, the indication should be checked again using ultrasound at the end of the drainage-free interval.

Example of catheter removal:

1. Check of the coagulation factors
2. Prepare a workstation under sterile conditions
3. Position the patient to ensure free access to the catheter
4. Generous skin infiltration with local anesthesia at the catheter exit site (observe exposure time)
5. Palpation of the polyester cuff from the outside (bulging thickening)
6. Incision of the skin, directly at the catheter outlet
7. Blunt, circular dissection of the polyester cuff using a button cannula, dissecting scissors, or similar (the catheter only grows together with the subcutaneous tissue* at the cuff; more difficult conditions may arise when removing the catheter due to the position itself or tumor walling).
8. Carefully place one hand on the tunnel section as an abutment, pull out the catheter with the other hand. If the catheter cannot be mobilized as a whole, an additional incision (corresponding to the length at the end of the tunnel section / the incision at the puncture site serves as a guide) is necessary in order to retrieve the catheter section lying/leading deep down there separately.
9. Final skin suture of the incision(s)
10. Apply a dressing.

GLOSSARY

Ascites – effusion or accumulation of fluid in the abdomen (stomach)

Protein loss – reduction in the amount of protein in the blood (protein depletion)

Electrolyte imbalance – disturbance of the balance of ions (e.g. sodium, potassium, calcium) in the blood or in the body cells

Empyema – infection of the pleura

Hemothorax – accumulation of blood in the pleural cavity

Hepatorenal syndrome – decrease in kidney function

Hypotension – low blood pressure

Infection – reaction of the body's own defenses to pathogens such as bacteria, viruses, etc

Circulatory collapse – temporary loss of consciousness

Peritoneal cavity – part of the abdominal cavity

Peritonitis – inflammation of the peritoneum (abdominal membrane)

Pleural effusion – effusion or accumulation of fluid in the pleural cavity

Pleurodesis – adhesion of the pleural cavity

Pneumothorax – air in the space between the lungs and the chest wall

Re-expansion pulmonary edema – fluid retention in the functional lung tissue due to the sudden expansion of the previously compressed lung; acute respiratory distress

Septa – tissue bridges / chambering between separated spaces in a body cavity (e.g. abdomen)

Subcutaneous tissue – tissue under the skin

Further, detailed information can also be found on our website at www.ewimed.de

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
EN This is the Summary of Safety and Clinical Performance (SSCP) location after the launch of European Database on Medical Devices/ EUDAMED

EN For the tearway sheath introducer is:

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EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague, The Netherlands

ewimed ewimed GmbH
Dorfstrasse 24
72379 Hechingen, Germany

ewimed Switzerland AG
Freiburgstrasse 453
3018 Bern
Switzerland

DE Behandlungseinheit gemäß Artikel 12, Absatz 3 der Richtlinie 93/42/EWG / **EN** Treatment unit in accordance with Article 12 (3) of Directive 93/42/EEC / **BG** Набор по смисъла на член 12, параграф 3 от Директива 93/42/ЕО / **DA** Behandlingscenter i overensstemmelse med artikel 12, stk. 3 i direktiv 93/42/EØF / **EL** Μονάδα θεραπείας σύμφωνα με το άρθρο 12 παράγραφος 3 της οδηγίας 93/42/ΕΟΚ / **ES** Unidad de tratamiento conforme al artículo 12, apartado 3 de la Directiva 93/42/CEE / **ET** Raviseade vastavalt direktiivi 93/42/EMÜ artikli 12 lõikele 3. / **FI** Hoitokeskus direktiivin 93/42/ETY 12 artiklan 3 kohdan mukaisesti / **FR** Unité de traitement conformément à l'article 12, paragraphe 3, de la directive 93/42/CEE / **HR** Komplet u skladu s člankom 12., stavkom 3. Direktive 93/42/EEZ / **HU** Kezelőközpont a 93/42/EGK iránylev 12. cikkének (3) bekezdésével összhangban / **IT** Kit procedurale in conformità all'articolo 12, paragrafo 3, della direttiva 93/42/CEE / **LT** Gydymo skyrius pagal Direktyvos 93/42/EEB 12 straipsnio 3 dalį / **LV** Procedūru komplekts saskaņā ar Direktīvas 93/42/EEK 12. panta 3. punktu / **MT** Unità ta' trattament skont l-Artikolu 12, il-paragrafu 3 tad-Direttiva 93/42/KEE / **NL** Behandelcentrum in overeenstemming met artikel 12, lid 3 van Richtlijn 93/42/EEG / **NO** Behandlingsenhet i henhold til artikkel 12, punkt 3 i direktiv 93/42/EØF / **PL** Ośrodek leczenia zgodnie z art. 12 ust. 3 dyrektywy 93/42/EWG / **PT** Centro de tratamento em conformidade com o n.º 3 do artigo 12.º da Diretiva 93/42/CEE / **RO** Unitate de tratament în conformitate cu articolul 12, punctul 3 din Directiva 93/42/CEE / **SK** Súprava pomôcok v súlade s článkom 12 ods. 3 smernice 93/42/EHS / **SL** Paket je skladen s členom 12(3) Direktive 93/42/EGS / **SV** Behandlingscenter i enlighet med artikel 12, punkt 3 i direktiv 93/42/EEG / **TR** 93/42/AEA sayılı direktifin 12. maddesinin 3. paragrafı uyarınca tedavi ünitesi

PleurX™/PeritX™ is a trademark and /or registered trademark of Becton Dickinson and Company, or one of its subsidiaries.

drainova™ is a trademark and /or registered trademark of ewimed.

octenisept® is a trademark of Schülke & Mayr GmbH.

ewimed

EW_10103Z_V11_2025-12

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EURASIA Deutschland GmbH
OBELIS S.A.
SISmed GmbH