
ewimed

ewimed Pneu-Pack II

DE / EN / BG / DA / EL / ES / ET / FI / FR / HR / HU / IT / LT / LV / MT / NL /
NO / PL / PT / RO / SK / SL / SV / TR

DE Gebrauchsanweisung / EN Instructions for use / BG Инструкции за употреба /
DA Brugsanvisning / EL Οδηγίες χρήσης / ES Instrucciones de uso /
ET Kasutusjuhend / FI Käyttöohjeet / FR Mode d'emploi / HR Upute za
uporabu / HU Használati utasítás / IT Istruzioni per l'uso / LT Naudojimo
instrukcijos / LV Lietošanas instrukcija / MT Struzzjoni għall-użu /
NL Gebruiksaanwijzing / NO Bruksanvisning / PL Instrukcja obsługi /
PT Instruções de utilização / RO Instrucţiuni de utilizare / SK Návod
na použitie / SL Navodila za uporabo / SV Instruktioner för
användning / TR Kullanım Kılavuzu

REF P7899

2 DE KONTAKT / EN CONTACT US / BG КОНТАКТ / DA KONTAKT / EL ΣΤΟΙΧΕΙΑ ΕΠΙΚΟΙΝΩΝΙΑΣ / ES CONTACTO / ET KONTAKTANDMED / FI YHTEYS / FR CONTACT / HR KONTAKT / HU KAPCSOLAT / IT CONTATTI / LT KONTAKTINĖ INFORMACIJA / LV KONTAKTINFORMĀCIJA / MT KUNTATT / NL CONTACT / NO KONTAKT / PL KONTAKT / PT CONTACTO / RO CONTACT / SK KONTAKT / SL KONTAKTNI PODATKI / SV KONTAKTA OSS / TR İLETİŞİM

AUSTRIA

ewimed austria GmbH
Tullnerbachstrasse 92a
3011 Neu-Purkersdorf
Tel.: +43 2231 2250-0
info@ewimed.com

BELGIUM

ewimed Belgium BV
Kortrijksesteenweg 1126a
9051 Sint-Denijs-Westrem
Tel.: +32 7848 0510
Info.be@ewimed.com

DENMARK

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ENGLISH

Read all the instructions for use carefully before using the product to drain an effusion.

All explanations of terms (*) can be found in the glossary on page 20.

THE EWIMED PNEU-PACK II

- ewimed Pneu-Pack II: for permanent drainage and drainage of effusions such as pleural effusion* and ascites*.
 - By gravity and for passive exhalation of air in pneumothorax* and pulmonary fistula*.
 - For patients from 10 kg body weight.
- Connection: Luer-Lock female
- Compatible irrigation tubes: drainova® irrigation tube (art. no. 6060, 6060/2) or PleurX™/PeritX™ irrigation tube (art. no. 50-7245A).
- Compatible catheters: drainova® / drainova® ArgentiC catheter (art. no. 1140 / 1150 / 1160 / 1240 / 1250 / 1260) or PleurX™/PeritX™ catheter (art. no. 50-7050 / 50-9050A).
- Minimum purchase 1 PU = 10 pieces, sterile
- Recommended additional product for drainage: drainova® dressing set (art. no. 8020).

The ewimed Pneu-Pack II is hereinafter referred to as the "reservoir".
The drainova® / drainova® ArgentiC catheter, PleurX™ / PeritX™ catheter is referred to below as the "catheter".

ADVANTAGES

Integrated flutter valve (Heimlich valve*) for maximum patient safety 1

- Prevents backflow of effusions and air into the body.

Air release valve 2

- Allows air to escape easily.

Effusion release valve 3

- Enables the effusion to be emptied cleanly and regularly.

Bed holder 4

- For simple and secure suspension of the reservoir.

WARNINGS

General

Do not use the product if:

- the sterile packaging is damaged
- the product is damaged
- the expiry date has passed

The reservoir is only suitable for use for the intended purpose.

The disposables are intended for single use only!

- The reservoir is supplied in a sterile condition, which is ensured by suitable quality controls.
- The reservoir must not be resterilised.
- The product is intended for single use only.
- Reuse may lead to cross-contamination.
- ewimed GmbH assumes no liability for products that have been reused, resterilised or removed from damaged packaging.

Drainage from the chest cavity

- Observe the drainage quantity according to the doctor's instructions.
- ewimed recommends draining no more than 1000 ml effusion per day.
- Only drain larger quantities according to the doctor's instructions.
- The patient may feel an unpleasant pulling sensation during drainage. If this is the case, close the pinch clamp on the connection tube to interrupt the flow of fluid.
- Possible complications of recurrent drainage of a pleural effusion* include pneumothorax*, re-expansion pulmonary oedema*, hypotension*, circulatory collapse*, septa* in the chest cavity or infection*.

Drainage from the abdominal cavity

- Observe the drainage quantity according to the doctor's instructions.
- ewimed recommends draining no more than 2000 ml effusion per day.
- Only drain larger quantities according to the doctor's instructions.
- The patient may feel an unpleasant pulling sensation during drainage. If this is the case, close the pinch clamp on the connection tube to interrupt the flow of fluid.
- Possible complications of recurrent drainage of ascites* include hypotension*, circulatory collapse*, electrolyte imbalance*, protein loss*, hepatorenal syndrome*, leakage of ascites* from incisions*, peritonitis*, septa* in the peritoneal cavity* and infection*.

Position of the reservoir during drainage 5

- Gravity-based drainage.
- The reservoir should always be positioned below the catheter exit site.
- Do not pull on the catheter exit site.
- Do not tug or pull on the connection tube.

Cleaning/disinfection

- In case of effusion leakage: Clean skin with soap and water.
- For surfaces: Use a suitable cleaning agent or disinfectant.
- Catheter / safety valve: Do not use octenisept® or disinfectant containing iodine.

The alcohol swabs are highly flammable.

Keep the wipes away from naked flames.

Change intervals of the reservoir

- The replacement interval is based on hygienic hospital standards.
- Change the reservoir every 72 hours at the latest.
- Check connection points 2x daily.
- In the event of accidental disconnection: Use a new sterile reservoir.

The reservoir is CE marked and sterilised with ethylene oxide gas (see packaging)

WARNINGS

In combination with the catheter

Never cut through the catheter or cut off the safety valve.

Avoid contact with sharp objects.

If the function of the safety valve is no longer guaranteed or if the safety valve has been accidentally cut off or is no longer present, proceed as follows: 6

- Press the catheter firmly closed with your fingers.
- Open the clamp with one hand and insert the catheter into the open end.
- Use a clamp available in the clinic. Position the clamp at a distance of approx. 2 cm from the valve.
- Contact the attending physician immediately.

Only use compatible products

- Only compatible products may be connected to the catheter.
- Do not insert anything other than the insertion pin specifically designed for this purpose into the safety valve, otherwise the safety valve may be damaged.
- Improper handling can lead to air entering the body or uncontrolled discharge of effusion.

Keep the safety valve on the catheter clean and the insertion pin of the irrigation tube on the connection tube sterile.

- Do not come into contact with non-sterile objects in order to avoid soiling or contamination.
- The safety valve must be disinfected before each drainage.

Make sure the safety valve and rinsing tube insertion pin are completely connected during the drainage process. Make sure the connection tube is not dragged or pulled.

- Ensure a sterile procedure when connecting the insertion pin to the safety valve.
- The plug connection must be firmly connected.
- Risk of contamination during disconnection. In such a case, disinfect the safety valve with a fresh alcohol swab.
- Dispose of disconnected drainage material.
- A new, sterile reservoir must be used for continuation.
- Observe the local disposal regulations.

Note: If required, the plug-in connection can be secured with a separately available drainova® clickFix (art. no. P1100 / P1100S) for additional protection against accidental disconnection.

The drainova® clickFix can be prescribed additionally and is not part of the reservoir. 7

CONNECTING THE PNEU-PACK II

- Connect the reservoir in accordance with the usual hygiene standards.

- Open the packaging and remove the reservoir.
 - Ensure that the Luer-Lock connection on the connection tube remains sterile.
- Remove a sterile drainova® irrigation tube or PleurX™ / PeritX™ irrigation tube from the packaging.
 - **drainova® irrigation tube:** Fully close the pinch clamp on the connection tube. 8
 - **PleurX™ / PeritX™ irrigation tube:** Close the roller clamp by turning it towards the male Luer lock connection. 9
- Connect the Luer-Lock female connection of the reservoir to the Luer-Lock of the irrigation tube. 10, 11
- Loosen the drainova® silicone cap, sterile or the PleurX™ / PeritX™ valve cap:
 - **drainova® silicone cap, sterile:**
 - Hold the safety valve of the catheter, carefully pull off the drainova® silicone cap, sterile. 12
 - Dispose of the cap.

• **PleurX™ / PeritX™ valve cap:**

- Hold the safety valve of the catheter, turn the PleurX™ / PeritX™ valve cap anti-clockwise and carefully it pull off. **13**
 - Dispose of the cap.
5. Hold the safety valve firmly with one hand.
- Clean the area around the safety valve opening (e.g. with an alcohol swab or a suitable alcohol-based disinfectant). **14**

Caution: Do not use octenisept® or disinfectants containing iodine.

• **drainova® irrigation tube:**

- Hold the insertion pin by the base.
 - Carefully remove the protective cover from the insertion pin with a slight twisting motion. **15**
 - Push the insertion pin into the safety valve.
 - Ensure that the safety valve and the insertion pin are firmly connected during the drainage process.
 - You should hear and feel them snap into place. **16**
- **PleurX™ / PeritX™ irrigation tube:**
- Hold the cap firmly and carefully pull the protective cover off the insertion pin. **17**
 - Push the insertion pin of the irrigation tube into the safety valve.
 - Secure the connection with a twisting motion.
 - You should hear and feel them snap into place. **18**
6. Open the pinch/roller clamp on the irrigation tube to start the drainage. **19, 20**
- Always position the reservoir below the catheter exit point, as drainage is by gravity. **5**
 - The effusion flows into the reservoir. The air release valve can be used for passive exhalation of air in pneumothorax* and pulmonary fistula*.
 - In the case of ascites* drainage, open or close the pinch/roller clamp at intervals in accordance with the doctor's instructions. **19, 20**

Note: The flow of fluid may be slowed or interrupted if:

- the effusion is almost completely drained
- the effusion has a high viscosity
- the safety valve is blocked by a coagulum* or a protein thread.

7. If necessary, open the air release valve.
- This valve is closed on delivery. **21**

Caution: If the reservoir is in a horizontal position, effusion may escape through the open air release valve.

8. Connect the bed holder to the reservoir.
- To hang up the reservoir, fold the bed holder backwards. **22**
 - Hang the holder on the bed frame, for example.
9. The permanent drainage is a safe and closed system.
- Change the reservoir/irrigation tube in accordance with the institution's hygiene regulations, but after 72 hours at the latest.
10. To empty the reservoir, open the effusion drain valve. **23**
11. When the desired drainage volume is reached, close the pinch/roller clamp. **19, 20**

• **PleurX™ / PeritX™ irrigation tube:**

- Hold the irrigation hose in one hand and the safety valve in the other.
 - Loosen the connection with a twisting motion and release the insertion pin of the irrigation tube from the safety valve with a firm, smooth motion. **25**
3. Clean the safety valve with an alcohol swab or suitable alcohol-based disinfectant. **14**
4. Attach the sterile drainova® silicone cap or the PleurX™ / PeritX™ valve cap:
- **drainova® silicone cap, sterile:**
 - Push the new sterile drainova® silicone cap firmly onto the safety valve of the catheter. **26**
 - The cap stays in position by itself without additional locking.
 - **PleurX™ / PeritX™ valve cap:**
 - Push the new valve cap onto the safety valve, turn it clockwise until it clicks firmly into place. **27**

GLOSSARY

Ascites – Effusion or accumulation of fluid in the abdomen (stomach)

Circulatory collapse – Temporary loss of consciousness

Electrolyte imbalance – Disturbance of the balance of ions (e.g. sodium, potassium, calcium) in the blood or in the body cells

Empyema – Infection of the pleura

Hepatorenal syndrome – Decrease in kidney function

Hypotension – Low blood pressure

Incision – Cut into the skin made during surgery, e.g. for implantation of a catheter

Infection – Reaction of the body's own defences to pathogens such as bacteria, viruses etc.

Integrated flutter valve (Heimlich valve) – Protection against backflow for increased patient safety

Peritoneal cavity – Part of the abdominal cavity

Peritonitis – Inflammation of the peritoneum

Pleural effusion – Effusion or accumulation of fluid in the pleural cavity

Pleurodesis – Adhesion of the pleural cavity

Pneumothorax – Air in the space between the lungs and the chest wall

Protein loss – Reduction in the amount of protein in the blood (protein depletion)

Pulmonary fistula – Open connection between lung tissue and pleural cavity

Re-expansion pulmonary oedema – Retention of fluid in the functional lung tissue due to the sudden expansion of the previously compressed lung; acute dyspnoea

Septa – Tissue bridges / chambering between separated spaces in a body cavity (e.g. abdomen)

Subcutaneous tissue – Tissue under the skin

You can also find more detailed information on our website at www.ewimed.com

FINAL STEPS

1. Open the packaging of the drainova® silicone cap, sterile or the PleurX™ / PeritX™ valve cap.
 2. Leave the cap in the packaging.
- **drainova® irrigation tube:**
 - Hold the connection tube in one hand and the safety valve in the other.
 - Pull the insertion pin out of the safety valve with a firm, smooth movement and put the connection tube down. **24**

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PleurX™/PeritX™ is a trademark and/or registered trademark of Becton Dickinson and Company, or one of its subsidiaries.
drainova® is a trademark and/or registered trademark of ewimed.
octenisept® is a trademark of Schülke & Mayr GmbH.
j. söllner AG
ewimed Switzerland AG
Cosanum AG