**REF 50-7500B, 50-7510, 50-7205B, 50-7210,  
50-0071, 50-7700, 50-9900C**

## **PleurX™ Drainage Procedure**

- en PleurX™ Drainage Procedure
- fr Procédure de drainage PleurX™
- de PleurX™ Drainagverfahren
- it Procedura di drenaggio con dispositivo PleurX™
- es Procedimiento de drenaje con PleurX™
- pt Procedimento de drenagem PleurX™
- nl PleurX™ drainageprocedure
- sv PleurX™ dräneringsprocedur
- fi PleurX™ dreneeraus
- no Prosedyre for PleurX™ drenasje
- da Procedure for dræning med PleurX™
- el Διαδικασία παροχέτευσης PleurX™
- tr PleurX™ Drenaj Prosedürü

**USA** Rx Only

For Single  
Use Only

Not made with  
natural rubber  
latex.

Do not use  
if package  
is damaged.

**DEHP**  
Contains phthalates

**STERILE** **EO**

**CareFusion**  
75 North Fairway Drive  
Vernon Hills, IL 60061 USA

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**Becton Dickinson Pty Ltd**  
Level 5, 66 Waterloo Road  
Macquarie Park, NSW 2113, Australia

**Becton Dickinson Limited**  
14B George Bourke Drive  
Mt. Wellington, Auckland, 1060 New Zealand  
362-24435 • 2020-12

© 2020 CareFusion Corporation or one of its affiliates.  
All rights reserved

For product inquiries: [bd.com/drainage](http://bd.com/drainage)  
PleurX, PeritX, CareFusion and the CareFusion logo  
are trademarks or registered trademarks of CareFusion  
Corporation, or one of its affiliates.

en

### Table of Contents

Product Information .....	1
How to Drain Using the PleurX™ Vacuum Bottle .....	2
Place a New Self-Adhesive Dressing .....	3
Dispose of the PleurX™ Drainage Bottle(s) .....	3
Frequently Asked Questions .....	4

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


























- Self-Adhesive Dressing
- (A) 
- (B) Blue Wrap around the following:  
Alcohol Pads
- (C) 
- (D) Gloves
- (E)   
Catheter Valve Cap
- (F) Blue Emergency Slide Clamp
- (F) 
- (G) Gauze Pads
- (G) 
- (H) Foam Catheter Pad
- (H) 

Flexible Bottle Cap  
Bouchon souple de la bouteille  
Flexibler Flaschendeckel  
Tappo flessibile del flacone  
Tapón flexible de la botella  
Tampa flexível do frasco  
Flexibele flesdop  
Bøjligt flasklock  
Taipuisa pullon kapseli  
Flexibel flaskehette  
Flexibelt flaskeålg  
Εύκαμπτο πώμα φιάλης  
Esnek Şişe Kapağı

Foil Seal  
Opercule  
Foliensiegel  
Sigillo laminato  
Sello de aluminio  
Vedante laminado  
Folieabdichtung  
Folieförlutning  
Lehtimetallikalvo  
Folieforsegling  
Φύλλο σφράγισης  
Folyo Kapak

Pointed end  
Extrémité pointue  
spitzes Ende  
estremità appuntita  
extremo en punta  
extremidade afiada  
puntig uiteinde  
spetsiga änden  
terävä pää  
spiss ende  
spids ende  
αχμηρό άκρο  
sivri ucunu

Vacuum Bottle  
Bouteille sous vide  
Vakuumflasche  
Flacone sottovuoto  
Botella de vacío  
Frasco de vacío  
Frasco de vacío  
Vacuümfles  
Vakuumflaska  
Tyhjiöpullo  
Vakuumflaske  
Φιάλη κενού  
Vakum Şişesi

Drainage Line  
Ligne de drainage  
Verbindungsschlauch  
Tubo di drenaggio  
Tubo de drenaje  
Linha de drenagem  
Drainagelijn

White "T" Plunger  
Piston blanc c en « T »  
Weißer T-Kolben  
Stantuffo bianco a "T"  
Êmbolo "T"  
Êmbolo em "T" branco  
Witte 'T'-zuiger  
Vit "T"-kolv  
Valkoinen T-mäntä  
Hvitt T-stempel  
Hvidt "T"-stempel  
Λευκό έμβολο σχήματος "T"  
Beyaz "T" Pistonu

Support Clip  
Pince de soutien  
Stützklamme  
Morsa di supporto  
Gancho de soporte  
Clipe de apoio  
Steunklem

Bottle Opener  
Dispositif d'ouverture du flacon  
Flaschenöffner  
Apriflacione  
Abridor de botellas  
Descapsulador de frascos


Rounded end  
Extrémité arrondie  
abgerundetes Ende  
estremità arrotondata  
extremo redondeado  
extremidade arredondada  
afgerond uiteinde  
rundade änden  
pyöreä pää  
avrundet ende  
afrundet ende  
στρογγυλεμένο άκρο  
yuvarlak ucunun

Access Tip Cover  
Capuchon de l'embout d'accès  
Hülle des Spezialansatzes  
Copertura punta di accesso  
Tapa de la punta de acceso  
Tampa da ponta de acesso  
Toegangstipdop  
Ätkomstspetsens hölje  
Liitinkärjen suojuus  
Deksel for tilgangsspiss  
Beskyttelse til adgangsspids  
Κάλυμμα άκρου πρόσβασης  
Giriş Ucu Kapağı

Access Tip  
Embout d'accès  
Spezialansatz  
Punta di accesso  
Punta de acceso  
Ponta de acesso  
Toegangstip

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## PleurX™ Drainage Procedure Instructions for Use


You should have been trained by a doctor or nurse in the correct way to use the items in this kit. If you have any problems or questions about draining the fluid, contact your doctor or nurse.

Read this entire pamphlet carefully before draining fluid.

### Product Information

The PleurX™ Catheter (**Figure 2**) may be placed in your chest to drain fluid from around your lungs or the PeritX™ Catheter may be placed in your abdomen to drain fluid that builds up in that area. Your doctor has placed the PleurX™ Catheter to allow you to drain fluid from your chest or the PeritX™ Catheter to allow you to drain fluid from your abdomen while you are at home. Part of the catheter is inside your body and has several holes that allow fluid to enter the catheter. The fluid is drained through the outer end of the catheter. There is a valve at the outer end of the catheter that prevents air from entering and fluid from leaking out when you are not draining.

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(2)

The PleurX™ Drainage Kit and Vacuum Bottles may be used to drain fluid through the PleurX™ Catheter or PeritX™ Catheter. To provide the best performance and to avoid damage to the catheter, always use the PleurX™ Drainage Supplies for draining fluid through the PleurX™ Catheter or PeritX™ Catheter.


### Indications for Use

The PleurX™ Drainage Kits are indicated for use only with the PleurX™ Catheter and PeritX™ Catheter for intermittent drainage.

### Warnings

**If you accidentally cut the catheter or the valve is leaking, follow these steps:**

- Pinch the catheter closed between your fingers.
- Slip the blue **emergency** slide clamp over the catheter and push the catheter completely into the small end of the clamp. This will close the catheter. (**Figure 3**)
- Notify your doctor immediately.

(3)

### Cautions

For single use only. Re-use may result in a non-functional product or contribute to cross contamination.

Contains Phthalates. The benefit of treatment outweighs the remote possibility of exposure to phthalates.

**Chest:** Potential complications of draining the pleural space include, but may not be limited to, pneumothorax, re-expansion pulmonary edema, hypotension, circulatory collapse and infection.

**Abdomen:** Potential complications of draining fluid from the abdomen include, but may not be limited to, hypotension, circulatory collapse, electrolyte imbalance, protein depletion, ascites leakage, peritonitis, wound infection, intraperitoneal adhesion and loculations of the peritoneal space.

### Contact your physician if:

- You believe your catheter is infected. Pain, redness (erythema), warmth to touch, swelling (edema), fever, or fluid around the catheter site may indicate your catheter is infected. Some discomfort and redness after insertion is expected but should not persist or worsen.
- Shortness of breath is not relieved after draining 1000 ml from your chest or 2000 ml of fluid from the abdomen at one time.
- You continue to experience symptoms, but little or no fluid drains from the catheter.
- Less than 50 ml of fluid drains in 3 drainage procedures in a row.
- The appearance (color, thickness, etc.) changes significantly between drainages.

### Sterility

The Procedure Pack and Vacuum Bottle (with drainage line) have been sterilized. The inside of the bottle and bottle connections are not exposed to sterilant. Samples should not be taken from the bottles. These products are for single use only. Do not re-sterilize.

Do not use if package is damaged.

CareFusion will not be responsible for any product that is re-sterilized, nor accept for credit or exchange any product that has been opened but not used.

### PleurX™ Drainage Supplies

The following supplies are required to complete a drainage procedure: (**Figure 1**)

- 1 Plastic Vacuum Bottle with attached Drainage Line
- 1 Procedure Pack containing the following sterile items:
  - 1 Self-Adhesive Dressing (**A**)
  - Blue Wrap (Wrapped around the following Items) (**B**)
    - 3 Alcohol Pads (**C**)
    - 1 Pair of Gloves (**D**)
    - 1 Valve Cap (**E**)
    - 1 Blue **Emergency** Slide Clamp (**F**)
    - 4 Gauze Pads 4" x 4" (10 cm x 10 cm) (**G**)
    - 1 Foam Catheter Pad (**H**)

## How to Drain Using the PleurX™ Vaccum Bottle

Now that your doctor or nurse has trained you on the drainage procedure, you should find it easy to follow these step-by-step instructions. You should drain as directed by your doctor, usually every one to two days. Do not change frequency or drain more fluid than your doctor has recommended without first consulting your doctor.

### Prepare to Drain

**Warning:** Do not use scissors or other sharp objects near the catheter.

**Note:** The blue **emergency** slide clamp should be used if you accidentally cut the catheter. **See Warnings (Figure 3)**

1. Set up a clean, clear workspace on a table or counter. Be sure to be near a sink and have a garbage can nearby.
2. Place the number of bottles that you plan on using near your workspace.
3. Thoroughly wash your hands with soap and water for at least 1 minute.
4. Remove the self-adhesive dressing from over your catheter and discard. Be sure not to tug on the catheter: support the skin while removing the self-adhesive dressing, gently grasp one edge and slowly peel the self-adhesive dressing from the skin. **(Figure 4a)** **Note:** Avoid skin trauma by peeling the self-adhesive dressing back, rather than pulling it from the skin. **(Figure 4b)**

**Caution:** Pain, redness (erythema), warmth to touch, swelling (edema), fever, or fluid around the catheter site may indicate your catheter is infected. Some discomfort and redness after insertion is expected but should not persist or worsen. If you see signs of infection, finish the drainage procedure and consult your doctor or nurse.

5. Thoroughly wash your hands again with soap and water for at least 1 minute.
6. Open the Drainage Kit bag and remove the Procedure Pack pouch.
7. Open the Procedure Pack pouch. Set the self-adhesive dressing aside.
8. Set packet with blue wrap on your workspace with the flap side up. Carefully unfold the blue wrap by pulling on the outside of the wrap. Leave the enclosed items on the wrap. The items and the inside of the wrap are sterile. Do not touch them with your ungloved hands or other non-sterile items. **(Figure 5)**
9. Remove the bottle with the drainage line from the outer pouch. Set the bottle near the blue wrap and place the access tip on the blue wrap. If the cover has fallen off of the access tip, make sure you do not touch the tip with your hands or anything non-sterile. It is especially important that the access tip remain sterile.
10. Remove the paper tape from the drainage line and uncoil the drainage line. Place the access tip back on the blue wrap.
11. Pick up the gloves by the folded cuff and put on as shown. **(Figure 6)** Both gloves fit either hand. Be careful not to let the outside of the gloves touch anything non-sterile, such as your skin or clothing.

12. Tear open the three alcohol pads, but do not remove the pads from their pouches. Place them on the blue wrap.

**Caution:** The alcohol pads are flammable. Do not expose the pads to an open flame.

### Connect the Drainage Bottle

**Warning:** Keep the valve on your catheter and the access tip on the drainage line clean. Keep them away from other objects to help avoid contamination.

**Warning:** Do not put anything except the access tip of the drainage line into the catheter valve, since any other device could damage the valve. A damaged valve may allow air into your body or let fluid leak out through the valve when you are not draining.

1. Close the roller clamp completely by rolling the wheel on the roller clamp toward the bottle. **(Figure 7)**

**Caution:** The roller clamp on the drainage line must be completely closed when not draining or the vacuum in the bottle may be lost.

2. Remove the cover from the access tip by twisting it and pulling gently. Discard the cover. Set the access tip back on the sterile blue wrap. **(Figure 8)**
3. Hold the base of the catheter valve and remove the cap by twisting it counterclockwise and pulling gently. **(Figure 9)** Discard the cap.
4. Clean around the valve opening with an alcohol pad. Discard the alcohol pad. **(Figure 10)**
5. Insert the access tip securely into the catheter valve. You will feel and hear a click when the access tip and valve are securely connected. **(Figure 11)**

**Caution:** Make sure that the valve and the access tip are securely connected when draining. If they are accidentally separated, they may become contaminated. If this occurs, clean the valve with an alcohol pad and use a new drainage line to avoid potential contamination.

**Caution:** Precautions should be taken to ensure the drainage line is not tugged or pulled.

### Drain Fluid

**Warnings:** It is normal to feel some discomfort or pain when draining fluid. If discomfort or pain is experienced when draining, roll the wheel on the roller clamp toward the bottle to slow or stop the flow of fluid for a few minutes. If you don't feel better after doing this, or the pain is severe, contact your doctor or nurse. Pain may be an indication of infection.

**Warnings:** Do not drain more than 1000 ml of fluid from your chest or 2000 ml of fluid from your abdomen at any one time.

1. The support clip on the vacuum bottle ensures that the vacuum seal on the bottle is not broken before you are ready. Remove the support clip by grasping the upper part of the flange and pulling outward in the direction of the arrow. **(Figure 12)**
2. Hold the bottle steady with one hand and push the white "T" plunger down with the other hand to puncture the foil seal. The vacuum in the bottle will pull the flexible bottle cap down. **(Figure 13)**

3. Roll the wheel on the roller clamp away from the vacuum bottle to begin drainage. **(Figure 14)** When fluid begins to flow into the drainage line you may partially close the roller clamp to slow the flow of fluid by rolling the wheel on the roller clamp toward the bottle.
4. The flow into the bottle may slow down when the fluid is almost completely drained. When the flow stops or the bottle is filled, completely close the roller clamp by rolling the wheel on the roller clamp toward the bottle. Drainage will usually take 5–15 minutes.

### Changing the Drainage Bottle

If you need to change to a new bottle for any reason, the procedure follows the same steps.

1. Remove the new bottle with the drainage line from the outer pouch. Set the bottle near the blue wrap and place the access tip on the blue wrap. If the cover has fallen off of the access tip, make sure you do not touch the tip with your hands or anything non-sterile. It is especially important that the access tip remain sterile.
2. Remove the paper tape from the drainage line and uncoil the drainage line. Place the access tip on the blue wrap.
3. Completely close the roller clamp by rolling the wheel on the roller clamp toward the bottle. **(Figure 7)**
4. Remove the cover from the access tip by twisting it and pulling gently. Discard the cover. Set the access tip back on the sterile blue wrap. **(Figure 8)**
5. To remove the used bottle, pull the access tip of the used bottle out of the valve in a firm, smooth motion. Set the used drainage line down.
6. Pick up the new drainage line and insert the access tip securely into the catheter valve. You will feel and hear a click when the access tip and valve are securely connected.
7. Resume following these instructions for use beginning with **Draining Fluid**. You do not need to clean the valve between bottles.

### Finish Drainage

1. Pull the access tip out of the valve in a firm, smooth motion. **(Figure 15)** Set the used drainage line down.
2. Clean the catheter valve with a new alcohol pad. Discard the alcohol pad.
3. Place the new cap over the catheter valve and twist it clockwise until it clicks into its locked position. Cap might feel loose until it is locked into place. **(Figure 16)**

### Place a New Self-Adhesive Dressing

1. Clean around the catheter site with a new alcohol pad. **(Figure 17)**
2. Place the foam catheter pad around the catheter. **(Figure 18)**
3. Wind the catheter into loops and place it over the foam pad. **(Figure 19)**

4. Cover the catheter with up to four (4) gauze pads. **(Figure 20)**
5. Remove gloves from both hands.
6. The self-adhesive dressing has three (3) layers:
  - a. printed liner
  - b. clear wound dressing
  - c. center panel and frame backing
7. Peel the printed liner from the self-adhesive dressing, exposing the adhesive surface. **(Figure 21)**
8. Center the self-adhesive dressing over the gauze pads and press it down. **(Figure 22)** **Note:** Do not stretch the self-adhesive dressing during application.
9. Remove and discard the center panel from the backing of the self-adhesive dressing. **(Figure 23)**
10. Slowly remove the frame while smoothing down the self-adhesive dressing edges. **(Figure 24)**
11. Smooth the entire self-adhesive dressing from the center toward the edges using firm pressure to enhance adhesion.

### Dispose of the PleurX™ Drainage Bottle(s)

1. Record the volume of drained fluid on drainage chart.

**Note:** Dispose of the drained fluid and the used drainage bottle per applicable local regulations. Used product may present a potential biohazard. Do not reuse.

If the applicable local regulations instruct you to empty the drainage bottle, follow steps 2-7 below.

2. Release any vacuum that may be left in the bottle by rolling the wheel on the roller clamp away from the bottle. Then completely close the roller clamp by rolling the wheel toward the bottle.
3. Place thumbs under the rounded end of the bottle opener and index fingers on top of the opposite end. Use thumbs to firmly push up to remove the bottle opener and flexible cap from the bottle. **(Figure 25)**
4. Remove the bottle opener from the drainage line by squeezing the flexible cap and pulling it through the bottle opener, or lifting the bottle opener over the drainage line. **(Figure 26)**
5. The bottle opener may be used to enlarge the hole in the foil seal to help fluid disposal. Hold the bottle steady with one hand and use your thumb to firmly press the pointed end of the bottle opener into the hole of the foil seal. Once inserted, rotate the bottle to make a larger opening in the foil seal. **(Figure 27)**
6. Remove the bottle opener from the bottle and empty the bottle into the toilet or sink.
7. Place the drainage line, bottle opener, and bottle in a plastic bag, seal tightly, and discard.

**Caution:** If fluid spills, use soap and water to clean skin and use appropriate cleaning agent for all other surfaces.

## Frequently Asked Questions

- Q: How often should I drain the fluid from my chest or abdomen?
- A: You should drain fluid as directed by your doctor, usually every one to two days. Consult your doctor before changing the frequency of your drainage.
- Q: What if I still feel short of breath or experience discomfort after I have finished draining?
- A: Notify your doctor if you continue to feel short of breath or experience discomfort. Do not drain more than 1000 ml of fluid from your chest or 2000 ml of fluid from your abdomen at any one time.
- Q: When will I know that the catheter can be removed?
- A: When you try to drain fluid 3 times in a row, and each time less than 50 ml drains into the bottle, you should see your doctor to find out if the catheter can be removed or if it needs to be replaced. (See next question.)
- Q: What does it mean if the volume of fluid I drain decreases or if I do not drain any fluid?
- A: If the fluid goes away suddenly or if the amount of drainage gradually declines, it is possible that the catheter or drainage line may be clogged. Squeeze the catheter and the drainage line gently. If drainage does not begin, follow the instructions for changing to another bottle. If the drainage does not start when you use a second bottle, call your doctor. If the amount of drainage gradually declines, the fluid may be drying up and it may be time for the catheter to be removed. Refer to question "When will I know that the catheter can be removed?"
- Q: How long will the catheter be in my chest or abdomen?
- A: The catheter will be in until fluid stops draining. The amount of time will vary from patient to patient. Fluid build-up is not likely to stop in the abdomen but may stop in the chest. The catheter may remain in place as long as you need it.
- Q: What should I do if the color of the fluid changes from the usual color?
- A: Any change in the appearance of the fluid should be reported to your doctor.
- Q: Can I take a shower or bath with the catheter in place?
- A: **Shower:** You can take a shower or sponge bath if a self-adhesive dressing like the one in the Procedure Pack is securely attached to your skin. The self-adhesive dressing is designed to keep fluid out. Be sure the self-adhesive dressing is completely and securely attached and that the catheter and gauze pads are all contained underneath it. If the gauze becomes wet when showering, remove the self-adhesive dressing immediately, clean and dry the area, and apply a new self-adhesive dressing as instructed under section "Place a New Self-Adhesive Dressing".
- A: **Bath:** Do not allow the catheter to soak underwater in a tub, bath, pool, etc.

Q: What happens if the PleurX™ Catheter or PeritX™ Catheter is accidentally pulled out?

A: In the unlikely event that the catheter is pulled out or the cuff becomes exposed, cover the exit site with a sterile self-adhesive dressing and seek immediate medical attention. The PleurX™ Catheter or PeritX™ Catheter has a polyester cuff that is normally under the skin where the catheter is inserted. The cuff and the sutures help keep the catheter in place.

Q: How will I know if my catheter is infected?

A: You should contact your doctor immediately if you believe your catheter is infected. Pain, redness (erythema), warmth to touch, swelling (edema), fever or fluid from around the catheter site may indicate your catheter is infected. Some discomfort and redness after insertion is expected but should not persist or worsen.

## Glossary of Terms

**Ascites:** liquid that builds up in the abdomen

**Pneumothorax:** air in the space between the lungs and the chest wall

**Circulatory collapse:** sudden drop in blood pressure: shock

**Erythema:** redness of the skin

**Edema:** swelling due to excessive fluid

**Electrolyte imbalance:** change in the normal level of certain chemicals in the body; such as sodium or potassium

**Hypotension:** low blood pressure

**Peritoneum:** membranes that line the inside of the abdomen and surround the abdominal organs

**Peritoneal space:** space in the abdomen formed by the inner and outer membranes of the peritoneum

**Peritonitis:** infection of the peritoneum

**Protein depletion:** reduction in the amount of protein in the blood

**Re-expansion pulmonary edema:** accumulation of fluid in the lung associated with very rapid drainage of a pleural effusion

## Accessories

Lockable Drainage Line      **REF 50-7245A**

Valve Cap      **REF 50-7235A**

Contact your doctor, home health care company, or CareFusion about how to obtain supplies.

## Warranty

CareFusion warrants that this medical device is free from defects in both materials and workmanship. The above warranties are in lieu of all other warranties, either expressed or implied, including any warranty of merchantability or fitness for a particular purpose. Suitability for use of the medical device for any surgical procedure shall be determined by a qualified medical professional. CareFusion shall not be liable for incidental or consequential damages of any kind.

**Note:** Not made with natural rubber latex.

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Becton Dickinson İth. İhr. Ltd.