
ewimed

DE drainova® Verbandmaterial-Set S / EN drainova® dressing material set S / BG Комплект превързочни материали S drainova® / DA drainova® forbindingsmateriale-sæt S / EL Σετ επιδεσμικού υλικού S drainova® / ES Juego de apósitos S drainova® / ET drainova® sidemete komplekt S / FI drainova® sidetarvikesetti S / FR drainova® kit de pansement S / HR Komplet zavojnog materijala S drainova® / HU drainova® kötszerkészlet S / IT drainova® kit di medicazione S / LT drainova® tvarstomosisios medžiagos rinkinys S / LV drainova® pārsienamo materiālu komplekts S / MT Sett ta' materjal ta' faxxatura S drainova® / NL drainova® verbandmateriaalset S / NO drainova® bandasjesett S / PL Zestaw materiałów opatrunkowych S drainova® / PT Conjunto de material de penso S drainova® / RO Set de pansamente S drainova® / SK Sada obvázového materiálu S drainova® / SL Komplet obvez S drainova® / SV drainova® uppsättning förbandsmaterial S / TR drainova® bandaj malzemesi seti S

DE Gebrauchsanweisung / EN Instructions for use / BG Инструкции за употреба / DA Brugsanvisning / EL Οδηγίες χρήσης / ES Instrucciones de uso / ET Kasutusjuhend / FI Käyttöohjeet / FR Mode d'emploi / HR Upute za uporabu / HU Használati utasítás / IT Istruzioni per l'uso / LT Naudojimo instrukcijos / LV Lietošanas instrukcija / MT Struzzjoni għall-użu / NL Gebruiksaanwijzing / NO Bruksanvisning / PL Instrukcja obsługi / PT Instruções de utilização / RO Instrucțiuni de utilizare / SK Návod na použitie / SL Navodila za uporabo / SV Instruktioner för användning / TR Kullanım Kılavuzu

REF 8020S

2 DE KONTAKT / EN CONTACT / BG КОНТАКТ / DA KONTAKT / EL ΣΤΟΙΧΕΙΑ ΕΠΙΚΟΙΝΩΝΙΑΣ / ES CONTACTO / ET KONTAKTANDMED / FI YHTEYS / FR CONTACT / HR KONTAKT / HU KAPCSOLAT / IT CONTATTI / LT KONTAKTINĖ INFORMACIJA / LV KONTAKTINFORMĀCIJA / MT KUNTATT / NL CONTACT / NO KONTAKT / PL KONTAKT / PT CONTACTO / RO CONTACT / SK KONTAKT / SL KONTAKTNI PODATKI / SV KONTAKTA OSS / TR İLETİŞİM

AUSTRIA

ewimed austria GmbH
Tullnerbachstrasse 92a
3011 Neu-Purkersdorf
Tel.: +43 2231 2250-0
info@ewimed.com

BELGIUM

ewimed Belgium BV
Kortrijksesteenweg 1126a
9051 Sint-Denijs-Westrem
Tel.: +32 7848 0510
Info.BE@ewimed.com

DENMARK

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ENGLISH

Read the entire instructions for use carefully before using the drainova® dressing material set S.

All explanations of terms (*) can be found in the glossary on page 16.

THE DRAINOVA® DRESSING MATERIAL SET S

- Intended use: The drainova® dressing material set S contains the necessary accessories for proper dressing of catheter exit sites and for absorbing body fluids after pleural effusion or ascites drainage.
- The drainova® silicone cap, sterile or PleurX™ / PeritX™ valve cap are used to seal these catheters.
- The set is sterilized with ethylene oxide gas.
- The sterile alcohol swabs are enclosed separately.

The drainova® dressing material set S is referred to below as the "dressing set".

We also recommend the following for changing the dressing set:

- Alcohol-based skin disinfectant spray
- Hand disinfectant
- Non-sterile disposable gloves

Indications:

The dressing material is intended for hygienic covering, fixation and padding of catheter exit sites and for absorbing body fluids.

Contraindications:

The dressing set has contraindications due to the components "foam slit dressing" and "polyurethane foam dressing", which can be found in the section "Precautions, warnings and recommended instructions for the foam slit dressing and the polyurethane foam dressing".

WARNINGS

General

The dressing set is only suitable for use for the intended purpose.

Target patient group:

The medical device may only be used on adults.

User group:

The medical device may only be used by adults.

Application environment:

The dressing material set is intended for use in clinical facilities (as well as hospitals, doctors' offices), care facilities and in the home environment.

Do not use the product if:

- the sterile packaging is damaged
- the product is damaged
- the expiry date of the set or the enclosed alcohol swabs has passed

The disposables are intended for single use only!

- The dressing set is delivered in a sterile condition, which is ensured by suitable quality controls.
- The dressing set must not be re-sterilized.
- Reuse may lead to cross-contamination.
- ewimed GmbH and the assembler Fuhrmann GmbH accept no liability for reused, re-sterilized products or products removed from damaged packaging.

Cleaning / disinfection

- In case of effusion leakage: Clean skin with soap and water.
- For surfaces: Use a suitable cleaning agent or disinfectant.
- Catheter / safety valve: Do not use octenisept® or disinfectant containing iodine.

Note: Any serious incident that occurs in connection with the use of this product should be reported to the manufacturer/assembler and the competent national authority of the Member State in which the user/patient is established.

Note: When disposing of used material, please observe the local, state and national regulations. If necessary, provide a residual waste container.

Precautions, warnings and recommended instructions for the foam slit dressing and the polyurethane foam dressing

The two products are summarised below as dressing material.

- These products are intended for single use and must not be reused. Repeated use can cause undesirable reactions such as itching and skin rashes and pose a risk of infection.
- If the packaging is damaged, the products must not be used.
- The dressing material should be used 2 cm beyond the edge of the wound.
- The dressing material should be checked regularly. If the wound secretion reaches within 2 cm of the absorbent edge, the dressing must be changed in good time.
- The frequency of dressing changes depends on the type of wound and the amount of wound fluid. (The dressing should be changed at least every 7 days; the total period of use must not exceed 30 days.)
- If the patient has a systemic or localised wound infection that is under control, however, these products may only be used under medical supervision.
- If the wound is caused by inadequate arterial circulation or a diabetic foot ulcer, these products may only be used under the guidance of a doctor or nurse and the wound must be checked daily.
- People who are allergic to important contact components of these products (e.g. silicone, polyurethane, acrylate, etc.) must not use them.
- These products must be removed before radiotherapy (including X-ray, ultrasound, hyperthermia and microwave therapy).

Please note:

- People who are allergic to polyurethane foam should use the products with caution.
- Not recommended for dry wounds.
- Use with caution in case of allergic disposition. In case of allergic reactions, discontinue use immediately and consult a doctor.

Note: It is strictly forbidden to use this product with hypochlorite solution, hydrogen peroxide and other oxidising agents, as these may affect the structure and performance of the product.

In combination with your catheter

Never cut through the catheter or cut off the safety valve.

Avoid contact with sharp objects.

If the function of the safety valve is no longer guaranteed or if the safety valve has been accidentally cut off or is no longer present, proceed as follows: **1**

- Press the catheter firmly closed with your fingers.
- Open the clamp with one hand and insert the catheter into the open end.
- Use a clamp available in the clinic.
Position the clamp at a distance of approx. 2 cm from the valve.
- Contact the attending physician immediately.

PREPARATION

17. Prepare a clean, clear work area.
18. Wash hands thoroughly with soap and water and then disinfect them.
19. If possible, put on non-sterile disposable gloves (not part of this set)
20. Remove the film dressing over the catheter.
- Carefully loosen one corner of the film dressing and slowly pull it off the skin.
- Make sure that the catheter is not accidentally pulled.

Note: In the event of swelling, signs of inflammation or fluid accumulation at the exit point of the catheter or at the removed skin interface, inform the attending physician.

21. If necessary, remove gloves and disinfect hands again.
22. Open the sterile packaging of the dressing set.

Note: You should provide an open container for the disposal of used material.

USING THE DRESSING SET

1. Place the dressing set on the work area with the flap facing upwards. **2**
 - Carefully unfold the surgical drape.
 - The items are sterile and must not come into contact with non-sterile objects.
2. Grasp the gloves at the opening and put them on. Do not touch any non-sterile objects! **3**
 - 2.1 Open the alcohol swab packaging, which is enclosed separately, on one side, and leave the swab in the packaging.
 - Place the swab on the edge of the surgical drape
3. Loosen the drainova® silicone cap, sterile or the PleurX™ / PeritX™ valve cap:
 - **drainova® silicone cap, sterile:**
 - Hold the safety valve of the catheter, carefully pull off the drainova® silicone cap, sterile.
 - Dispose of the cap. **4**
 - **for PleurX™ / PeritX™ valve cap:**
 - Hold the safety valve of the catheter, turn PleurX™ / PeritX™ valve cap counterclockwise and carefully pull off.
 - Dispose of the cap. **5**
4. Clean the safety valve with an alcohol swab. **6**
 - Observe exposure time of 30 seconds.

Caution: Do not use octenisept® or disinfectants containing iodine.

5. Open the packaging of the drainova® silicone cap, sterile and drop it onto the surgical drape. **7**
6. Attach a new drainova® silicone cap, sterile and push it firmly onto the safety valve of the catheter. **8**
 - The cap holds independently without additional locking.
7. Clean the area around the catheter exit point with the alcohol-based skin disinfectant spray. **9**
 - Allow a contact time of 30 seconds.
 - Do not apply the new dressing until completely dry.
8. Place the foam slit dressing around the catheter. **10**
9. Remove the protective film from the fixation strip on the side facing the catheter and fix the foam slit dressing to the skin. **10**
10. Remove the protective film from the upper side (side facing away from the skin) and attach the silicone cap or protective cap to the adhesive strip. **10**
11. Remove gloves.
12. Remove the two coated strips on the adhesive side (side facing the skin). **11**
13. Hold the film dressing at the edge strips.
14. Centre the film dressing by placing the PUR foam over the foam slit dressing. **12**
15. Remove the two-part cover film (on the side facing away from the skin). **13**
16. Finally, smooth out the film dressing again.

Note: Dispose of the used dressing material in accordance with local regulations.

16 FREQUENTLY ASKED QUESTIONS

When can the stitches be removed?

The skin stitch can be removed after 7 - 10 days. However, the retaining suture should only be removed after 30 days.

Can I shower or bathe with the catheter?

If the catheter has been completely and correctly covered with the self-adhesive film dressing, you can take a shower. If the foam slit dressings get wet while showering, remove the dressing, clean and dry the area and apply a new dressing. The catheter must not be submerged in water in a bath or swimming pool. Avoid direct contact of the valve with bath or shower water.

Note: We recommend showering before the drainage / dressing change so that the dressing can be changed promptly if it gets wet.

How do I find out if I have an infection*?

Inform a doctor immediately if you experience pain, reddening of the skin (erythema*), tangible warmth, swelling (oedema*), fever, or fluid accumulation, as these may be signs of an infection*.

Proceed in the same way if changes in the colour of the effusion become visible.

Redness and slight pain are to be expected shortly after implantation, but these should not persist or worsen.

What should I do if I accidentally pull out the catheter?

The catheter should fuse with the subcutaneous* skin tissue through the polyester sleeve. If the catheter is accidentally pulled out and the cuff is

visible, cover the catheter exit site with a sterile dressing and contact the implanting clinic immediately.

RETURN OR CREDIT NOTE

Only originally packed (unopened, unlabelled, and undamaged) boxes of dressing sets will be accepted for return and credited.

GLOSSARY

Ascites – Effusion or accumulation of fluid in the abdomen (stomach)

Erythema – Reddening of the skin

Infection – Reaction of the body's own defences to pathogens such as bacteria, viruses, etc.

Oedema – Swelling due to excessive fluid accumulation

Pleura – Thin, serous skin in the chest cavity (pleura)

Subcutaneous tissue – Tissue under the skin

For more detailed information visit our website at www.ewimed.de

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PleurX™/PeritX™ is a trademark and/or registered trademark of Becton Dickinson and Company, or one of its subsidiaries.
drainova® is a trademark and/or registered trademark of ewimed.
octenisept® is a trademark of Schülke & Mayr GmbH.
Albo-Healthcare GmbH
ewimed Switzerland AG