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Patient leaflet about the implant

This information leaflet contains important instructions for the safe use of your implant. Please keep this document together with your implantation card and present it during medical examinations.

1. The drainova® catheter

Intended purpose: The catheter is used as an implant product for recurrent, refractory, malignant and non-malignant effusion accumulations in serous body cavities to alleviate symptoms by draining the fluids. This includes: malignant pleural effusion, non-malignant pleural effusion, malignant ascites, non-malignant ascites.

This means: The catheter is implanted when fluid repeatedly accumulates in a patient's body and other treatments are no longer sufficiently effective. It is used to drain this fluid from specific body cavities - particularly the chest cavity (around the lungs) or the abdominal cavity - in order to relieve associated discomfort. It is used for both cancer-related and non-malignant conditions.

The implanted catheter is MRI-safe. MRI examinations may be performed under the specified conditions. Please inform medical staff about your implant before every MRI examination.

2. Information on the service life of the implant

The product is intended for a duration of use of more than 30 days depending on the patient's clinical condition. Therefore, the duration of implantation of the drainova® catheter is not limited.

3. Follow-up care and check-up

Your treating physician will determine the frequency and type of appointments based on your individual needs.

Next follow-up appointment Physician's stamp

Date

4. Material and components

The implant consists of the following essential materials, which may be relevant for allergies or intolerances:

Silicone tube: Silicone (vinyl dimethyl polysiloxane), radiopaque material barium sulfate (BaSO_4)

Cuff: Polyester fleece (polyester), silicone adhesive (dimethyl polysiloxane with deblocked acetoxo groups)

Safety valve: Acrylonitrile butadiene styrene (ABS), polypropylene (PP), polyisoprene (PI), polydimethylsiloxane

Please inform your healthcare professional if you have any known allergies or material intolerances (e.g., to metals, silicone, plastics, latex, etc.).

5. Instructions for safe everyday use

Please follow these recommendations for safe use:

Always carry your implant card with you.

- + Inform healthcare professionals about your implant.
- + Follow your physician's specific instructions.
- + Only the recommended compatible products designated for drainage should be used.
- + Do not pull on the catheter.
- + Observe the maximum recommended drainage volume.
- + Keep the safety valve on the catheter clean and the insertion pin on the connection tube sterile. Make sure the safety valve and insertion pin are fully connected during the drainage process. Make sure the connection tube is not dragged or pulled.
- + In case of effusion leakage: Clean skin with soap and water. For surfaces: Use a suitable cleaning agent or disinfectant. Catheter/safety valve: Do not use octenisept® or disinfectant containing iodine.
- + The catheter must not be submerged in water in a bath or swimming pool. Avoid direct contact of the safety valve with bath or shower water.
- + The duration of implantation of the catheter is not limited. Clinical data from post-market surveillance for an implantation duration of up to 30 days and further data from post-market surveillance for an implantation duration of up to 270 days are available. In individual cases, the implantation time may be shortened due to the clinical situation and a medical assessment (e.g. infection).

6. Contact

If you have any questions, complaints, or problems, please contact your healthcare professional or manufacturer at any time. Please refer to the contact details on page 40.

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